

KIDDIES OMEGA POWER



NAME	OMEGAA1 POWER
PRODUCT CODE	A1POWK
SEAT	14"
CASTORS	8X2" SOLID
REAR TYRE	12½ X 2¼ " FOAM FILLED
CONTROL	DYNAMIC REM D01
BATTERY	12V36A/H
ARMREST AND FOOT	REMOVABLE & HEIGHT ADJUSTABLE ARM REMOVABLE FOOTREST

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

ALPHA A6 POWER WITH TILT



NAME	ALPHA TILT POWER
PRODUCT CODE	A6POW
SIZE	12", 14"
BATTERIES	2 X 36AH/12 VOLT
MOTORS	250 WATT
CASTORS	8 X 2" FOAM FILLED GREY
SPEED	6KM/H
RANGE	20KM (TERRAIN AND USER DEPENDENT)
REAR RIMS & TYRES	121/2 X 21/4 ALUMINIUM RIM 121/2 X 21/4 ALUMINIUM RIM
ARMS	HEIGHT ADJUSTABLE T ARM
GRADIENT	12 DEGREES
CONTROLS	DYNAMIC SHARK REMD11 POWER MODULE PMB11
SEAT FUNCTION	
LOAD CAPACITY	LOAD CAPACITY – 120KG WEIGHT – 70KG
PACKING	70X70X100CM

DRAGON JUNIOR POWER WHEELCHAIR

The *Invacare Dragon Junior* is a powerful and adjustable power wheelchair to meet everyday needs of the child. The multiple seat adjustments and the motor options make this wheelchair easy to use.



Technical datas

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Seat angle	Backrest angle
<i>Dragon Junior</i>	290-360 mm 320-400 mm	300 - 360 mm	425 - 475 mm	450 mm	180 - 250 mm	200 - 350 mm	0° à 17° (powered or man.)	From -10° to 30° (man.) or from 5 to 30° (powered)
	Driving unit width	Total height	Total length inc. footrests	Total length without footrests	Total weight	Max. user weight	Battery capacity	Turning diameter
<i>Dragon Junior</i>	610 mm	950 mm	1030 mm	850 mm	100 kg	90 kg	2 x 60 Ah	1660 mm
	Max. Climbable obstacle height	Max. safe slope	Motor capacity	Driving range	Speed			
<i>Dragon Junior</i>	60 mm	18 %	2 x 180 W	30 km	6 km/h			

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za