

TRANSPORTER BUGGY



Weight – 12kg
Load Capacity – 50kg
Sizes – 12", 14", 16"
Frame – Mild Steel
Rear Tyre – 11" Solid

NAME	BUDDY OFFROAD BUGGY
PRODUCT CODE	BUG001
SIZE	12", 14", 16"
FRAME	MILD STEEL
TYRE	11" SOLID REAR WHEEL
LOAD CAPACITY	50KG
WEIGHT	12KG

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

BUDDY BUGGY URBAN



Item Code Yk8003

Frame Powder coating
Front wheel 8"×1.25" ,PU tires
Rear wheel 12" PU tires
Handgrip PVC
Seat and backrest Nylon and woven cloth, Depth and width fixed
Front fork Aluminum
Footplate Plastic
Anti-tipper tube
Quick release for rear wheel
Footrest height adjustable
Breathable headrest and backrest
High-quality sunshade"
Seat Width 14"

Carton size: 1300*420*360mm
GW/NW: 18.5/16.0KG

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INT ROADSTER BUGGY



Yk8001 which seat width is 280mm-320mm
 YK8002 which seat width is 320mm-400mm

- Aluminium frame, light weight and luxury design
- Foldable chair frame for easy transportation and storage
- Comfortable soft seat, back over and adjustable soft headrest
- Height adjustable & foldable Pushbar
- Removable solid base seat system
- Seat Back height adjustable and seat base depth adjustable
- Reclining Backrest
- Adjustable angle & height footrest
- Mesh bag under for small belongings

Seat width 280mm-320mm YK8001
 Seat width 320mm-400mm YK8002
 Color: Orange & grey
 Blue & grey



FOLDED & SEAT REMOVED

FOLDED:	Specification 82cm*66cm*45cm	Yk8001	20kg
FOLDED:	Specification 82cm*74cm*45cm	Yk8002	22kg

CPC WITH TILT AND RECLINE



Model ALK958LC-36

- 1) Aluminium lightweight frame with comfortable seat, reclining high back, adjustable headrest
- 2) 6" PVC castor 16" wheel
- 3) Adjustable and detachable armrest, elevating footrest, seat and backrest adjustable
- 4) With elasticity belt

NAME	CPC WITH TILT AND RECLINE
PRODUCT CODE	ALK958LC / YK9140
SIZE	12", 14", 16", 18"
FRAME	ALUMINIUM LIGHTWEIGHT
TYRE	6" PVC CASTOR 16" WHEEL
LOAD CAPACITY	50KG (18" - 60KG)
WEIGHT	15KG

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE AZALEA MINOR



PRODUCT CODE: BNC0014

Rea Azalea Minor – designed for petite adults and teenagers
 The Rea Azalea Minor is designed for users that require a seat that has a smaller depth and width. Adapted from the Rea Azalea, the Azalea Minor boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism. This model is ideal for all users that require a solid, stable wheelchair. Now with new frame colours and upholstery. A tension adjustable backrest is available as a new option.



Tension adjustable backrest
 A new option adjustable in height and as wide as the seat.



Table tray
 It is transparent and the rail offers gripping opportunities for posture solutions.



Shorter Chair
 The chassis is 30 mm shorter, the seat depth is 50 mm shorter and the backrest 100 mm lower than on standard models. Here in new colour, electric green.

Technical data

Rea Azalea Minor	340 - 440 mm	380 - 450 mm	400/450 mm	550 - 350 mm	240 - 360 mm	330 - 500 mm	+1° to 25°
Rea Azalea Minor	0° - 30°	SW + 220 mm ISW + 230 mm	900 - 1250 mm	900 - 970 mm	32 kg	Max 75 kg	Transport weight 14,5 kg
	<p>Rea Azalea Minor complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.</p>						<p>Invacare International Sàrl Route de Cité-Ouest 2 - 1196 Gland Switzerland Tel: +41 (0)22 354 60 10 Fax: +41 (0)22 354 60 11 hqeuropa@invacare.com www.invacare.eu.com</p>
Rea Azalea Minor	Crash-tested ISO 7176-19						

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za

FREEDOM DESIGN NXT



Transfers as easy as A, B, C...

- A:** Travels light, weighing just 22 pounds.
- B:** Folds into small package.
- C:** Fits into small spaces.

Seat Frame Widths

The NXT is available in 11 sizes from 10" wide to 20" wide.
 Each size from 10" wide to 16" wide adjusts an additional 2" in width.
 The 17" wide size adjusts an additional 1" in width to 18" wide.
 The 18", 19" and 20" wide sizes are set.

Seat Frame Depths

The NXT is available in 3 sizes: a short, a medium or a long.
NEW NXT Mini: For x-short 8" to 12" range.
 Each size adjusts an additional 4" in depth.
 Seat frame depth is measured without seating, from the front of the back post to the front end of the seat frame.
 The short range is 12" to 16."
 The medium range is 14" to 18."
 The long range is 16" to 20."

Seat Depth Extensions

For additional seat frame depths, extensions can be added to the back of the seat frame. There is a 2", a 3" or a 4" extension to choose from.

Rider Capacity

225 lbs.

Frame

Available in 25 colors.

Seat to Floor Height

The NXT can be configured in a variety of seat to floor heights.
 Selections in front casters and rear wheel sizes determine a range of choices from 12.5" low to 18" high, and most everything in between.

Tilt Range

45 degrees.

TIAN FEEDING & POSTURE CHAIRS

PRODUCT CODE: TC401 & TC403



Code: TC401

SPECIFICATION

Dimensions(L*W*H)mm:630*430*(840~950).

Seat size(L*W*H)mm:330*(300~360)*360.

Load capacity: ≥ 70 kg

Net weight: 16kg

Adopting carbon dioxide welding, electrostatic spraying, anti-aging and good anti-rust.

All pad are a foam forming.

Distance can be adjust of pillow height: 0-140mm

Gear pad adjustment: 0-100mm

Pedal height adjustment:0-80mm.

Plastic moulded tray, quick & simple dismount for easy of use

The distance of front-back wheels: diameter ≥ 50 mm, Adjustable.

SPECIFICATION

- Table to seat cushion distance 210mm
- Head pad lift adjustment range 0~200mm
- Crotch pad before and after the adjustment range 0~100mm
- Seat to the maximum distance of the pedal 200mm
- Load Capacity 40kg
- Overall dimensions (L×W×H)680×560×800



Code: TC403

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za