

BASIK STANDING

PRODUCT CODE ALK158



Model /**ALK158**
Name /*Electric Wheelchair*



- 1) Spray powder frame with breathable mesh cushion
- 2) 8" big castor, 12" PU rear wheel
- 3) Detachable plate, the armrest can be turned backwards
- 4) can go up and down and walk upright by electric motor
- 5) with safe light
- 6) Seat width of 40cm

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

ROC HYDRAULIC STAND-UP



Innovate, simple, light weight design to stand-up and sit-down

- * Advanced stability design of safety concern
- * Unique easy control make stand-up and sit-down possible
- * Dual safety protection to prevent accidental improper operation
- * Multiple adjustments to fit user size

NAME	ROC (SMC-S0) HYDRAULIC STAND-UP
PRODUCT CODE	49311
CONTROLS	DUAL GAS ASSISTED LIFT PUMPS
CASTORS	6X1" GREY CASTORS
REAR RIM	24X1 3/8" PNEUMATIC GREY TYRE ON ALUMINIUM RIM
BATTERIES	
LOAD CAPACITY	LOAD CAPACITY - 110KG WEIGHT - 30KG
ADJUSTABILITY	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY
PACKING SIZES	65X68X94CM 14", 16", 18"

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

GRIFFIN SEMI-ELECTRIC STAND-UP



NAME	SEMI-ELECTRIC STAND-UP
PRODUCT CODE	50003
CONTROLS	SINGLE ACTUATOR, FINGER OPERATED
CASTORS	6X1" GREY CASTORS
REAR RIM	24X1 3/8" PU GREY TYRE ON ALUMINIUM SPOKE RIM
BATTERIES	2 X SEALED 3.8AH/24VOLT
LOAD CAPACITY	LOAD CAPACITY - 110KG WEIGHT - 30 KG
ADJUSTABILITY	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY
PACKING	65X68X94CM
SIZES	14", 16", 18"

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE DRAGON VERTIC

Access all areas – with the push on a button
The Invacare® Dragon Vertic is an agile power chair that not only provides mobility to an active user, but will bring you face to face with other people. Its comfortable seating system easily moves you into a standing position by a simple joystick command, providing increased access and excellent mobility.



NAME	INVACARE DRAGON VERTIC
PRODUCT CODE	VERTIC
CONTROLS	SHARK II STEERING
CASTORS	MINI CATORS
REAR RIM	PUNCTURE PROOF REAR TYRES
STD ACCESSORIES	
BATTERIES	2X12V 60aH
LOAD CAPACITY	LOAD CAPACITY – 131KG WEIGHT – 100KG
ADJUSTABILITY	ARM REST, HEAD REST, SLIDING BACK

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

PHOENIX STAND-UP POWER



NAME	SATURN PHOENIX POWER
PRODUCT CODE	49312
SIZE	14", 16", 18"
BATTERIES	2XSEALED 36AH/12 VOLT
CASTORS	8X2"FOAM FILLED GREY TYRE
REAR TYRES	3.00-8 FOAM FILLED GREY TYRE
CONTROLS	DYNAMIC SHARK REMD11 DUAL FUNCTION SHARK REMOTE REMD11 POWER MODULE PMB11
LOAD CAPACITY	LOAD CAPACITY – 110KG GROSS WEIGHT – 90KG
ADJUSTABILITY	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT
STD ACCESSORIES	SOLID FIXED CLOTHES GUARDS PERSPEX
PACKING	103X72X72CM

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

EASYADAPT EVOLV SIT-TO STAND FRAME

Modular design makes the EasyStand Evolv the most versatile standing frame available. The base unit is functional for many users, however when more support is needed a wide variety of positioning options can be added to accommodate the most involved users.

Medical Benefits of Standing

- Improve/maintain range of motion
- Decrease joint/muscle contractures
- Management of atrophy in the trunk and leg muscles
- Improve strength to trunk and lower extremities
- Decrease muscle spasms
- Improve/maintain bone integrity/skeletal development
- Lessen/manage the progression of scoliosis
- Manage pressure (ulcers) through changing positions
- Improve bowel function and regularity
- Aid in kidney and bladder functions
- Strengthen cardiovascular system and build endurance
- Improve circulation
- Reduce swelling in lower extremities



Standard Front
with options



Shadow Tray
with options

NAME	EASY STAND EVOLV SIT-TO-STAND
PRODUCT CODE	PNG50084
CONTROLS	MEDICAL HAND PUMP
CASTORS	6" SOLID
STD ACCESSORIES	SEAT – BACK AND TRAY UPCHARGE OPTIONS
LOAD CAPACITY	LOAD CAPACITY - 115KG
ADJUSTABILITY	SEAT DEPTH

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

TIAN PARALLEL BARS

PRODUCT CODE: KR-PXG-03



SPECIFICATIONS:

Length 3,000mm
Height: 1,120mm (Adjustable 780-1,220) Width: Adjustable

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za