

# INVACARE HEAVY DUTY TRANSPORT CHAIR



## NAME

INVACARE  
HEAVY DUTY  
TRANSPORT  
CHAIR

## PRODUCT CODE

TRHD22FR

## SEAT WIDTH

22"

## SEAT DEPTH

18"

## BACK HEIGHT

20"

## SEAT-TO-FLOOR HEIGHT

19.5"

## OVERALL WIDTH OPEN

28"

## OVERALL DEPTH W/O RIGGINGS

36.5"

## LOAD CAPACITY

400LB

## Features

- 13 month limited warranty on side frames and cross braces
- 12" rear wheels with flat free tyres
- Removable padded desk arms
- Carbon steel frame with dual cross braces
- Swingaway footrest with aluminium footplates
- Flame retardant black nylon upholstery

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# SATURN HERCULES – XL



<b>NAME</b>	SATURN HERCULES
<b>PRODUCT CODE</b>	HER022
<b>SIZE</b>	22"
<b>COLOURS</b>	COBOLT BLUE
<b>FRAME</b>	REINFORCED POWDER COATED MILD STEEL
<b>ARMS</b>	REMOVABLE DESK OR FULL ARM
<b>FOOTREST</b>	REMOVABLE, ELEVATING LEG REST OPTIONAL EXTRA
<b>CASTORS</b>	8x2" PU CASTOR ON 5MM STEEL FORK
<b>RIMS</b>	24X1 <sup>3/8</sup> MAG RIM WITH 24X1.75 PU TYRE - FIXED AXLE
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON
<b>STRENGTH</b>	SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL
<b>LOAD CAPACITY/ PACKING</b>	LOAD CAPACITY 125KG WEIGHT - 21KG 100 X 35 X 85 (CM)

# JADE BARIATRIC



**Specifications:**  
PRODUCT CODE: JADE

Size	22", 24", 26"
Seat Depth	18"
Frame	Reinforced Powder Coated Mild Steel
Footrest	Removable
Upholstery	Black Reinforced Nylon
Load Capacity	145kg – 185kg
Optional	Elevating leg rest

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# INVACARE BARIATRIC IVC 9000 TOPAZ



HDPG option 1, 000 Fb  
Product weight (9TPZA)  
Limited warranty 3 years

## **Specifications:**

PRODUCT CODE: 9TPZ

### Seat Sizes:

Widths 20", 22", 24", 26" 28", 30"

Depths 18", 20"

Seat-to-floor heights 17.5", 19.5"

15.5" with different wheel and caster packages

Back heights – adjustable 16"-18"

Front riggings

Swingaway footrest

Swingaway padded elevating legrest

### Arm styles

Full-length fixed height, padded, removable

Desk-length fixed height, padded, removable

Full-length adjustable height, padded, removable

Desk-length adjustable height, padded, removable

(All arm styles available conventional)

Weight capacity 700 lb

HDPG option 1,000 lb

Product weight (9TPZ/WD2818/28) 85lb

Limited warranty 3 years on frame and cross braces

The wheelchair is designed specifically to suit the unique needs of bariatric patients. The Topaz combines the strength and durability of a triple chrome-plated, carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty cross braces, this chair provides not only ease in manoeuvrability but stability that will be trusted for years to come.

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# INVACARE ACTION 4

## Where stability meets style!

A comfortable, foldable, lightweight wheelchair, the *Invacare® Action<sup>4</sup><sub>NG</sub>* is ideally suited for everyday use.

A key feature of this next-generation model is the innovative Invacare® hallmark of cross-compatibility; a shared platform with other models in the Action family as well as other wheelchairs in our range means that users and advisors can easily update or refurbish their wheelchairs and enjoy maximum opportunities for modularity and personalized adjustment.



### Options

**NEW**



**Heavy Duty Version 160 kg**  
A wider choice of seat widths (55.5 & 60.5 cm), reinforced chassis and a new greater carrying capacity (maximum user weight 160 kg).




**Hemi version\***  
A 38.5 cm seat height, dual handrim propulsion, hemi armrest and built out castor housing.



**Amputee Legrest**  
The amputee legrest is both height- and angle- adjustable.

### Technical data




38, 40.5, 43, 45.5, 48, 50.5, 55.5, 60.5 cm



40, 45 and 50 cm



38.5\*, 41, 43.5, 46, 48.5, 51 cm



38 - 46 cm  
39 - 51 cm\*\*



20 - 23 - 26 cm,  
24 cm, 19 - 29.5 cm




35 - 48 cm



0° - 6°



7°  
0° - 30°



Seat width + 21 cm




90 - 103 cm



105 - 111 cm



83 - 90 cm



30.5 cm



15.2 kg



125 / 160 kg



Transport:  
Weight 8.3 kg



Crash tested  
ISO 7176-19\*\*\*

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# MATRIX VI CUSHION



## Invacare® Matrix® Vi Cushion

Model no. MA(WWDD)-Vi

### Features and Functional Benefits

#### Outer cover:

- Moisture resistant and breathable
- Includes slip resistant base, hook and loop fasteners, lifting strap
- Available in black and Kid\*ab\*ra print

#### Inner cover:

- Moisture resistant
- Easily wipes clean
- Zippered cover allows access for cushion modifications

#### Top layer:

- Soft (ILD) High Resilient (HR) foam for optimal comfort and support
- Pre-ischial contour provides stability to pelvis
- Abductor/adductor contour promotes midline leg positioning
- Made with Ultra-Fresh™ for antibacterial and odor protection

#### Visco foam insert:

- Molded visco-elastic foam ischial relief insert
- Anatomically targeted to ischial/sacral area
- Provides immersion to sensitive bony prominences

#### Bottom layer:

- Firm (ILD) High Resilient (HR) foam for support and durability
- Trochanteric shelf provides stability to hips
- Waterfall front edge
- Rear radius for sacral support

**HCPCS Code**  
**E2607/E2608**

### STANDARD

Width	14", 16", 18", 20", 22"
Depth	15" – 20"
A. Height at adductor	3.75"
B. Height at abductor	3.25"
C. Height at leg trough	2.5"
Weight capacity	300 lb.



Kid\*ab\*ra fabric standard on pediatric sizes

### HEAVY-DUTY

Width	18" – 30"
Depth	18" – 26"
A. Height at adductor	3.75"
B. Height at abductor	3.25"
C. Height at leg trough	2.5"
Weight capacity	600 lb.

### PEDIATRIC

Width	12", 14"
Depth	12", 14"
A. Height at adductor	2.875"
B. Height at abductor	2.75"
C. Height at leg trough	2"
Kid*ab*ra fabric standard	
black fabric optional	

### DROP BASE



Includes built-in drop base, drop hooks, DuraClip fasteners  
Allows adjustment of seat angle

### PELVIC OBLIQUITY KIT



Easily installed on all Invacare Matrix Vi/ Kid\*ab\*ra-Vi cushions/ bases  
Accommodates .5" or 1" obliquities

### CUSTOM MODIFICATIONS



Leg length discrepancies, width/ depth reductions  
Foam, inner and outer covers modified to your specifications

# BARIATRIC H/DUTY ELECTRIC BED



Function:	Electric Hi-Lo 3 Function
Castors:	5" PA Castors (foot w/brake)
Frame:	Reinforced Steel – Powder Coated
Outer Dimensions:	2225x1070mm
Surface Dimension:	2020x880mm
Adjustable Height:	460-960mm
Bed Head Adj:	To 70 degrees G.W: 136kg
Packing:	Foot: 117x23x112cm (53kg)
	Head: 112x23x110cm (42kg)
	H&F Board: 112x19x70 (34kg)
Upper Side Rail	74x7x45cm (7kg)
Includes:	Trunk Side Rail Set, Emergency Manual Crank
Load Capacity:	247kg

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za