

S542 YOGA FOLDING SCOOTER

PORTABLE
CONVENIENCE


TRAVEL COMPANION




Easy-to-read
control panel



Small folding dimension
for easy storage



Patented Twin-front-wheels with
spring suspension system



Optional suitcase

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

CTM HS268 FOLDING SCOOTER MANUAL



TECHNICAL SPECIFICATIONS

Overall Length	953 mm / 37.5"
Overall Width	500 mm / 19.7"
Overall Height	960 mm / 37.8"
Front Wheels	200 mm / 8"
Rear Wheels	200 mm / 8"
Weight W/ Batteries	28.6 kg / 63.1 lbs
Weight Of Heaviest Piece	22.7 kg / 50 lbs
Max. Speed	6.4 kmph / 4 mph
Weight Capacity	115 kg / 255 lbs
Ground Clearance	43 mm / 1.7"
Grade Climbable	6 degree
Curb Climbable	40 mm / 1.6"
Turning Radius	1120 mm / 44.1"
Brake	Electro-Mechanical
Seat Type	Padded Foldable
Seat Width	419 mm / 16"
Motor Size	250W, 3800 r.p.m
Battery Size	(1) 24V. 11.6Ah
Battery Weight	2.8 kg / 6.2 lbs
Travel Range	20 km / 12.4 Miles
Battery Charger	2A Off Board
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Button

Features



CTM HS268 FOLDING SCOOTER AUTO



TECHNICAL SPECIFICATIONS

Overall Length	953 mm / 37.5"
Overall Width	500 mm / 19.7"
Overall Height	960 mm / 37.8"
Front Wheels	200 mm / 8"
Rear Wheels	200 mm / 8"
Weight W/ Batteries	29.5 kg / 65 lbs
Weight Of Heaviest Piece	23.6 kg / 52 lbs
Max. Speed	6.4 kmph / 4 mph
Weight Capacity	115 kg / 255 lbs
Ground Clearance	43 mm / 1.7"
Grade Climbable	6 degree
Curb Climbable	40 mm / 1.6"
Turning Radius	1120 mm / 44.1"
Brake	Electro-Mechanical
Seat Type	Padded Foldable
Seat Width	419 mm / 16"
Motor Size	250W, 3800 r.p.m
Battery Size	(1) 24V, 11.6Ah
Battery Weight	2.8 kg / 6.2 lbs
Travel Range	20 km / 12.4 Miles
Battery Charger	2A Off Board
Electronics	Sensor Key, Touch Screen, Coded Lock, Remote Control for Auto-Folding

Features



Medop^{cc}



Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

