

CTM HS115 MICRO SCOOTER



HS-115			
Type	Mobility Scooters	Curb Climbing	35 mm / 1.4"
Model Number	HS-115	Turning Radius	1090 mm / 43"
Overall Length	1010 mm / 39.8"	Supension	N/A
Overall Width	520 mm / 20.5"	Brake	Electro-Mechanical
Overall Height	905 mm / 35.6"	Seat Type	Swivel Padded Foldable
Wheels: Front	195 mm / 8"	Seat Width	425 mm / 17"
Wheels: Rear	195 mm / 8"	Motor Size	250W, 4700 r.p.m.
Weight w/ Batteries	37 kg / 82 lbs	Battery Size	(2) 12V, 12Ah
Max. Speed	6.4 kmph / 4 mph	Battery Weight	9.2 kg / 20 lbs
Weight Capacity	115 kg / 255 lbs	Travel Range	10 km / 6.2 Miles
Ground Clearance	30 mm / 1.2"	Battery Charger	2A Off Board
Grade Climbable	6 degree	Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob
Color	Red		

* Subject to change without notice.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za

CTM HS118 MICRO SCOOTER



Type	Mobility Scooters	Curb Climbing	35mm/1.4"
Model Number	HS-118	Turning Radius	1270mm/50"
Overall Length	1000mm/39.4"	Suspension	N/A
Overall Width	510mm/20"	Brake	Electro-Mechanical
Overall Height	890mm/35"	Seat Type	Padded Foldable
Wheels: Front	195mm/8"	Seat Width	425mm/17"
Wheels: Rear	195mm/8"	Motor Size	250W, 4700 rpm
Weight w/ Batteries	40kg/88lbs	Battery Size	(2) 12V • 12Ah
Max. Speed	6.4kmph/4mph	Battery Weight	9.2kg/20lbs
Weight Capacity	115kg/253lbs	Travel Range	10km/6.2miles
Ground Clearance	30mm/1.2"	Battery Charger	1.8A Off Board
Grade Climbable	8 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Blue		

CTM HS265 SCOOTER



HS-265			
Type	Mobility Scooters	Curb Climbing	50 mm / 2"
Model Number	HS-265	Turning Radius	1060 mm / 41.7"
Overall Length	1020 mm / 40.2"	Suspension	N/A
Overall Width	500 mm / 20"	Brake	Electro-Mechanical
Overall Height	890 mm / 35"	Seat Type	Swivel Padded Foldable
Wheels: Front	200 mm / 8"	Seat Width	425 mm / 16.7"
Wheels: Rear	200 mm / 8"	Motor Size	200W, 4650 r.p.m
Weight w/ Batteries	40.1 kg / 88.4 lbs	Battery Size	(2) 12V, 12Ah
Max. Speed	8 kmph / 5 mph	Battery Weight	9.5 kg / 20.9 lbs
Weight Capacity	136 kg / 300 lbs	Travel Range	11 km / 6.8 Miles
Ground Clearance	40 mm / 1.6"	Battery Charger	2A On Board
Grade Climbable	8 degree	Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob
Color	Burgundy, Blue		

* Subject to change without notice.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za

CTM HS295 AND HS295U SCOOTER



FEATURES:

- Maximum speed: 8kmph
- Capacity: 136kg
- Wireless design for trouble free dismantling
- Combination battery and charger pack
- Standard Delta handlebar

TECHNICAL SPECIFICATIONS

Overall Length	1010 mm / 39.8"
Overall Width	550 mm / 21.7"
Overall Height	900 mm / 35.4"
Front Wheels	210 mm / 8 1/4 "
Rear Wheels	210 mm / 8 1/4 "
Weight W/ Batteries	50.3 kg / 110.8 lbs
Max. Speed	8 kmph / 5 mph
Weight Capacity	130 kg / 287 lbs
Ground Clearance	60 mm / 2.4"
Grade Climbable	8 degree
Curb Climbable	40 mm / 1.6"
Turning Radius	1300 mm / 51.2"
Brake	Electro-Mechanical
Seat Type	Swivel Padded Foldable
Seat Width	445 mm / 17.5"
Motor Size	200W, 4650 r.p.m
Battery Size	(2) 12V. 18Ah
Battery Weight	13.4 kg / 29.5 lbs
Travel Range	16 km / 9.9 Miles
Battery Charger	2A Off Board
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob



FEATURES:

- Maximum speed: 8kmph
- Capacity: 130kg
- Wireless design for trouble free dismantling
- Combination battery and charger pack
- Two 12V.18Ah batteries for longer travel range
- Standard Delta handlebar

CTM HS328 - MID SIZE SCOOTER



Large Seat Option



- Maximum speed: 8kmph
- Capacity: 136kg
- Easy take apart mechanism
- Patented battery pack locking mechanism
- Standard LED headlight
- 22Ah battery pack for longer travel range
- Standard Delta handle bar
- Optional tiller head charging socket

Medop cc
















Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE COLIBRI


PRODUCT CODE: 1578998



Technical data

				
Invacare Colibri 465 mm	400 mm	370 - 420 mm	385 mm	505 mm
				
Invacare Colibri 1010 mm	44 kg (12 Ah batteries) 49 kg (18 Ah batteries)	136 kg	2 x 12 Ah 2 X 18 Ah	2200 mm
				
Invacare Colibri 45 mm	6° / 10.5%	200 W	Up to 16 km	8 km/h

Shroud colours



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

SCOOTER ACCESSORIES



Canopy



Sun Shade



Oxygen Tank Holder



Cane Holder



22" Captains Seat



Flat Free Tyre

Medop^{cc}



Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

