

CTM HS320 MID-LINE SCOOTER



Large Seat Option



Solid Tyre Option



CTM MID-LINE SCOOTER – HS320

| Type | Mobility Scooters | Curb Climbing | 50mm/2" |
|----------------------------|--------------------------------------|------------------------|---|
| Model Number | HS-320 | Turning Radius | 1020mm/40.2" |
| Overall Length | 1090mm/42.9" | Suspension | N/A |
| Overall Width | 550mm/21.7" | Brake | Electro-Mechanical |
| Overall Height | 870mm/34.3" | Seat Type | Light Weight Foldable Swivel |
| Wheels: Front | 230mm/9" | Seat Width | 460mm/18" |
| Wheels: Rear | 230mm/9" | Motor Size | 450W 3700 rpm |
| Weight w/ Batteries | 53.8kg/118.7lbs | Battery Size | (2) 12V • 18 Ah |
| Max. Speed | 8kmph/5mph | Battery Weight | 12.7kg/28lbs |
| Weight Capacity | 135kg/300lbs | Travel Range | 135km/8.4 miles |
| Ground Clearance | 40mm/1.6" | Battery Charger | 2A Off Board |
| Grade Climbable | 12 degrees | Electronics | On/Off Key Switch, Battery Level Indicator, Speed Control Knob |
| Colour | Burgundy, Silver, Blue, Turquoise | | |

CTM HS360 MID-LINE SCOOTER



CTM MID-LINE SCOOTER – HS360

| Type | Mobility Scooters | Curb Climbing | 50mm/2" |
|---------------------|--------------------------------------|-----------------|---|
| Model Number | HS-360 | Turning Radius | 1415mm/55.7" |
| Overall Length | 1160mm/45.7" | Suspension | N/A |
| Overall Width | 550mm/21.7" | Brake | Electro-Mechanical |
| Overall Height | 870mm/34.3" | Seat Type | Light Weight Foldable Swivel |
| Wheels: Front | 230mm/9" | Seat Width | 460mm/18" |
| Wheels: Rear | 230mm/9" | Motor Size | 450W 3700 rpm |
| Weight w/ Batteries | 67kg/148.9lbs | Battery Size | (2) 12V • 18Ah |
| Max. Speed | 8kmph/5mph | Battery Weight | 12.7kg/28lbs |
| Weight Capacity | 135kg/300lbs | Travel Range | 13km/8.1 miles |
| Ground Clearance | 40mm/1.6" | Battery Charger | 2A Off Board |
| Grade Climbable | 12 degrees | Electronics | On/Off Key Switch, Battery Level Indicator, Speed Control Knob |
| Colour | Burgundy, Silver, Blue, Turquoise | | |

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

CTM HS570 GT-3 SCOOTER



Solid Tyre Option

CTM GT-3 SCOOTER – HS570

| Type | Mobility Scooters | Curb Climbing | 50mm/2" |
|---------------------|-----------------------------------|-----------------|--|
| Model Number | HS-570 | Turning Radius | 1080mm/42.5" |
| Overall Length | 1240mm/48.8" | Suspension | N/A |
| Overall Width | 610mm/24" | Brake | Electro-Mechanical |
| Overall Height | 965mm/38" | Seat Type | Premium Upholstery Mid-Back w/ Headrest Swivel Seat-Sliding Adjustment |
| Wheels: Front | 260mm/10" | Seat Width | 455mm/18" |
| Wheels: Rear | 260mm/10" | Motor Size | 500W 3400 rpm |
| Weight w/ Batteries | 84kg/185.2lbs | Battery Size | (2) 12V• 36 Ah |
| Max. Speed | 8kmph/5mph | Battery Weight | 24.7kg/55lbs |
| Weight Capacity | 135kg/300lbs | Travel Range | 30km/18.6 miles |
| Ground Clearance | 40mm/1.6" | Battery Charger | 5A Off Board |
| Grade Climbable | 12 degrees | Electronics | On/Off Switch, Battery Level Indicator, Speed Control Knob |
| Colour | Burgundy, Silver, Blue, Turquoise | | |

CTM HS580 SUPREME SCOOTER



CTM SUPREME SCOOTER – HS580

Solid Tyre Option

| Type | Mobility Scooters | Curb Climbing | 50mm/2" |
|----------------------------|-----------------------------------|------------------------|--|
| Model Number | HS-580 | Turning Radius | 1520mm/59.8" |
| Overall Length | 1280mm/50.4" | Suspension | N/A |
| Overall Width | 610mm/24" | Brake | Electro-Mechanical |
| Overall Height | 960mm/37.8" | Seat Type | Premium Upholstery Mid-Back w/ Headrest Swivel Seat-Sliding Adjustment |
| Wheels: Front | 260mm/10" | Seat Width | 455mm/18" |
| Wheels: Rear | 260mm/10" | Motor Size | 500W 3400 rpm |
| Weight w/ Batteries | 87.9kg/194lbs | Battery Size | (2) 12V• 36 Ah |
| Max. Speed | 8kmph/5mph | Battery Weight | 24.7kg/55lbs |
| Weight Capacity | 135kg/300lbs | Travel Range | 30km/18.6 miles |
| Ground Clearance | 40mm/1.6" | Battery Charger | 5A Off Board |
| Grade Climbable | 12 degrees | Electronics | On/Off Switch, Battery Level Indicator, Speed Control Knob |
| Colour | Burgundy, Silver, Blue, Turquoise | | |

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za

CTM HS585 SCOOTER



Solid Tyre Option

| HS-585 | | | |
|----------------------------|-----------------------------------|------------------------|--|
| Type | Mobility Scooters | Curb Climbing | 50 mm / 2" |
| Model Number | HS-585 | Turning Radius | 1480 mm / 58.3" |
| Overall Length | 1280 mm / 50.4" | Supension | N/A |
| Overall Width | 610 mm / 24" | Brake | Electro-Mechanical |
| Overall Height | 1110 mm / 43.7" | Seat Type | Swivel Mid-Back W/ Headrest & Seat Sliding Mechanism |
| Wheels: Front | 260 mm / 10" | Seat Width | 455 mm / 18" |
| Wheels: Rear | 260 mm / 10" | Motor Size | 500W, 3400 r.p.m |
| Weight w/ Batteries | 89.3 kg / 196.9 lbs | Battery Size | (2) 12V. 36Ah |
| Max. Speed | 8 kmph / 5 mph | Battery Weight | 24.7 kg / 55 lbs |
| Weight Capacity | 135 kg / 300 lbs | Travel Range | 34 km / 21 Miles |
| Ground Clearance | 40 mm / 1.6" | Battery Charger | 5A Off Board |
| Grade Climbable | 12 degree | Electronics | On/Off Key Switch, Battery Level Indicator, Speed Control Knob |
| Color | Burgundy, Silver, Blue, Turquoise | | |

* Subject to change without notice.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE LEO

PRODUCT CODE: 1579613



Features and Accessories



Finite tiller movement
Users can easily adjust the tiller to suit their needs.



Brake light
Rear brake lighting warns other people that the scooter is slowing down, thus preventing collisions.



Comfort seating
Swivel functionality and sliding rails ensure optimal seating position.



Two-step disengaging lever
This unique concept prevents the scooter from accidentally freewheeling.

Technical data



470 mm



410 mm



460 - 560 mm



475 mm



590 mm



1220 mm



59 kg (excl. batteries)
83 kg (incl. batteries)



136 kg



2 x 36 Ah



2620 mm



60 mm



70 mm



10°



240 W



38 km



8 km/h

* Seat in maximum fifth rearward locking position from the front.

Colour of shrouds



Moonstone Silver



Onyx Blue



Ruby red

Upholstery



Black vinyl

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

CTM BARIATRIC SCOOTER



PRODUCT CODE: H5890



Solid Tyre Option

| Type | Mobility Scooters | Curb Climbing | 70mm/2.8" |
|---------------------|-----------------------------------|-----------------|--|
| Model Number | HS-890 | Turning Radius | 1730mm/68.1 |
| Overall Length | 1530mm/60.2" | Suspension | Front & Rear |
| Overall Width | 685mm/27" | Brake | Electro-Mechanical |
| Overall Height | 1070mm/42.1" | Seat Type | Sliding and Swivel Reclining |
| Wheels: Front | 320mm/13" | Seat Width | 480mm/19" |
| Wheels: Rear | 320mm/13" | Motor Size | 700W 5400 rpm |
| Weight w/ Batteries | 114kg/251lbs | Battery Size | (2) 12V • 50Ah |
| Max. Speed | 15kmph/9.3mph | Battery Weight | 31.3kg/69lbs |
| Weight Capacity | 227kg/500lbs | Travel Range | 37.5km/23.4 miles |
| Ground Clearance | 85mm/3.3" | Battery Charger | 5A Off Board |
| Grade Climbable | 10 degrees | Electronics | On/Off Key Switch, Battery Level Indicator, Speed Control Knob |
| Colour | Burgundy, Silver, Blue, Turquoise | | |

HITCH RACK X-RAD TRAILER



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za