

WHEELCHAIR ACCESSORIES AND OPTIONS

SINGLE PIECE RAMPS



PRODUCT CODE: RAM

Aluminium Frame
 Sizes - 36", 48", 72"
 Single piece foldable with hinge and handle
 250kg Load Capacity
 Width - 34"

TELESCOPIC RAMPS



PRODUCT CODE: RAM020

15cm Internal width
 25cm External width
 3 Section Aluminium Frame
 Telescopic 110cm to 280cm
 Gross Weight per set - 25kg
 5cm Lip Height
 Rough Grip along the surface and on the lip
 Load Capacity - 200kg

CHIN CONTROL WITH SEAT



TRANSFER BOARD



PRODUCT CODE: TRB002

EZ ENABLER TRAY



PRODUCT CODE: 10858

ABS table for active wheelchairs
 Aluminium clamp with e-clip x 2
 With bag
 Packing: Box size: 43 x 25.2 x 13cm

WHEELCHAIR ACCESSORIES AND OPTIONS

**PORTABLE O2 HOLDER
(ALSO A MODEL WITH DRIPSTAND)**



WHEELCHAIR TRAY



TAB (TENSION ADJUSTABLE BACK)

“T”: DRIP ROD



**ELEVATING LEG REST FOR
STEEL WHEELCHAIR
(HEEL STRAP SOLD SEPARATELY
ON SELECTED WHEELCHAIRS ONLY)**

DESK AND FULL ARM



WHEELCHAIR ACCESSORIES AND OPTIONS

ANTI-TIP WHEEL



24 X 1" SHOX TYRE & 5" CASTOR



8 X 1" SOLID CASTOR



**AMPUTEE BRACKET: 5CM
CENTER OF GRAVITY ADJUSTMENT**



8 X 2" PU OR FOAM FILLED CASTOR



**24X2" REAR RIM WITH
12.7MM BEARING**



**5X1 1/4" & 5X1" CASTOR
ALSO AVAILABLE IN 6"**



WHEELCHAIR ACCESSORIES AND OPTIONS

HEAD REST WITH HARDWARE



4 POINT HARNESS



**LAP/BODY STRAP
(ALSO WITH VELCRO)**



**MATRX ELITE
DEEP BACK**



**MATRX EASY SET
HARDWARE**



**BODYPOINT PIVOT
SHOULDER HARNESS**



**BODYPOINT
PELVIC BELT**

