

SATURN ALU-LITE



NAME	SATURN ALU-LITE
PRODUCT CODE	LW25-77
SIZE	14", 16:, 18" AND 20"
COLOURS	METALLIC BLUE, WET BLACK
FRAME	FOLDING ALUMINIUM POWDER COATED
ARMS	FIXED HEIGHT FLIP BACK DESK ARM
FOOTREST	REMOVABLE
CASTORS	8X1" SOLID CASTOR ON ALUMINIUM FORK
RIMS	24X1 ^{3/8} MAG RIM WITH PU TYRE - 12MM QR AXLE
UPHOLSTERY	BLACK REINFORCED NYLON TAB
ADJUSTMENT	CENTRE OF GRAVITY, LOW SLUNG, HEIGHT ADJUSTABLE BACKREST WITH PUSH HANDLES
OPTIONS	ELEVATING LEG REST, HEIGHT ADJUSTABLE ARMS, 8X2" CASTOR ON STEEL FORK, ALUMINIUM SPOKE RIM, TENSION ADJUSTABLE BACK UPHOLSTERY, AMPUTEE BRACKET, ANTI-TIP WHEELS
WEIGHT	LOAD CAPACITY 115KG GROSS WEIGHT 14 KG
PACKING	100 X 34 X 85 (CM)

ROC HYDRAULIC STAND-UP



Innovate, simple, light weight design to stand-up and sit-down

- * Advanced stability design of safety concern
- * Unique easy control make stand-up and sit-down possible
- * Dual safety protection to prevent accidental improper operation
- * Multiple adjustments to fit user size

NAME

ROC (SMC-S0) HYDRAULIC STAND-UP

PRODUCT CODE

49311

CONTROLS

DUAL GAS ASSISTED LIFT PUMPS

CASTORS

6X1" GREY CASTORS

REAR RIM

24X1 3/8" PNEUMATIC GREY TYRE ON ALUMINIUM RIM

LOAD CAPACITY

LOAD CAPACITY - 110KG

WEIGHT - 30KG

ADJUSTABILITY

BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY

PACKING

65X68X94CM

SIZES

14", 16", 18"

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

GRIFFIN SEMI-ELECTRIC STAND-UP



NAME	SEMI-ELECTRIC STAND-UP
PRODUCT CODE	50003
CONTROLS	SINGLE ACTUATOR, FINGER OPERATED
CASTORS	6X1" GREY CASTORS
REAR RIM	24X1 3/8" PU GREY TYRE ON ALUMINIUM SPOKE RIM
BATTERIES	2 X SEALED 3.8AH/24VOLT
LOAD CAPACITY	LOAD CAPACITY - 110KG WEIGHT - 30 KG
ADJUSTABILITY	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY
PACKING	65X68X94CM
SIZES	14", 16", 18"

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za



Invacare Action Series

INVACARE Action Family: The changing face of Mobility

The Action Family product offer consists of:

Action2 - Standard aluminium wheelchair
Action3 Junior - Paediatric wheelchair from age 3 and up

Action3 - Lightweight wheelchair, single crossbar, highly configurable with a huge range of options, including:

- Hemiplegic version
- One arm drive version
- Double hand rim version
- Comfort version
- and more

Action4 - Lightweight wheelchair, double crossbar, highly configurable with huge range of options, including:

- Hemiplegic version
- One arm drive version
- Double hand rim version
- Comfort version

Heavy-duty option up to 160kg user weight



Cross Compatibility table: Action² Action³

Leg rest	x	x
Footplates	x	x
Armrest	x	x
Arm pads	x	x
Side guards	x	x
Seat upholstery	x	x
Backrest canes	x	x
Backrest upholstery		x
Push handles	x	x
Brakes	x	x
Fork		x
Castor wheels	x	x
Frame	x	x
User weight	125KG	125KG

Action³ One arm drive version



Action³ Comfort version



Action³ Hemiplegic version



INVACARE ACTION 2

The Next Generation has arrived!

The all-new **Action²_{NG}** has been developed from the tried-and-tested **Action³** concept. The specifications meet the basic needs of the users whilst refusing to compromise on safety or comfort. The next generation **Action²_{NG}** shares the same essential platform as the **Action³**, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the **Action²_{NG}** perfectly fits into the Action family.



Features



Swing in and out footrest
The 80 degree footrest swings in and out for easy transfers in or out of the chair.



Adjustable center of gravity
The newly-designed rear wheel plate can be adjusted by 5 cm.



Anti-tippers
Provide extra safety on slopes or uneven ground.

Technical data

38, 40,5, 43, 45,5, 48 cm	40,5 / 43 cm	48,5, 51 cm	43 cm	20 - 26 cm	35 - 48 cm	80° / 0 - 80°	7° / 0° - 30°
SW+ 20,5 cm	93,5 cm	82 / 84,5 cm	41 cm	15 kg	max. 125 kg	Crash tested ISO 7176-19*	Crash tested ISO 7176-19*

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE ACTION 3

A star is re-born!

The **Action³NG** is one of the best examples of cross-compatibility in the Invacare® range. Offering high levels of efficiency for dealers, platform sharing across the Action family and other Invacare® models makes re-use and refurbishment an easy and attractive feature for users. An appealing 'next generation' design together with improved functionalities in terms stability and comfort means that the **Action³NG** marks the rebirth of our flagship product in style.



Features



14,2 kg weight
The total weight of 14,5 kg addresses the need for improved handling, storage and transportation. The lowest transportation weight (without legrest, armrest and rear wheels) is 7,6 kg.



Brake and Brake extension
Fold the brake handle for easier side transfers.



Length compensating Leg support
New elevating legrest offers optimal length compensation with the elevation stepped in approx. 3 – 5 degrees.



Footrest
The proven design of the swing in / swing out footrest has an easy-to-operate quick release mechanism. Adjust the foot support height in steps of 1 cm.



20 rear wheel positions
Driving performance is enhanced with a choice of 20 rear wheel positions as standard. Users can set the rear wheel to Standard, Medium, Active or Amputee position.



Magnesium castor attachment
A castor attachment provides ease of angle adjustment and enhanced support. The use of the magnesium attachment reduces the weight of the chair by almost 400 gram.

Technical data




38, 40,5, 43, 45,5, 48, 50,5 cm




40, 45 cm



38,5, 41, 43,5, 46, 48,5, 51 cm



40, 43, 51 cm



20 - 23 - 26 cm, 24 cm, 19 - 29,5 cm




35 - 48 cm, 34 - 47 cm



80° / 0° - 80°



0° / 7° / 0° - 30° / -6° - +24°



Seat width + 20,5 cm




81,5 - 102 cm




100 cm



72 - 85 cm



28,5 cm



14,2 kg



max. 125 kg



Transport weight is 7,6 kg (without armrest, footrest and rearwheels)



Crash tested ISO 7176-19*

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE ACTION 4

Where stability meets style!

A comfortable, foldable, lightweight wheelchair; the *Invacare*® **Action**⁴_{NG} is ideally suited for everyday use.

A key feature of this next-generation model is the innovative *Invacare*® hallmark of cross-compatibility; a shared platform with other models in the Action family as well as other wheelchairs in our range means that users and advisors can easily update or refurbish their wheelchairs and enjoy maximum opportunities for modularity and personalized adjustment.



Options

NEW



Heavy Duty Version 160 kg

A wider choice of seat widths (55.5 & 60.5 cm), reinforced chassis and a new greater carrying capacity (maximum user weight 160 kg).



Hemi version*


A 38.5 cm seat height, dual handrim propulsion, hemi armrest and built out castor housing.



Amputee Legrest

The amputee legrest is both height- and angle- adjustable.

Technical data




38, 40.5, 43, 45.5, 48, 50.5, 55.5, 60.5 cm




40, 45 and 50 cm



38.5*, 41, 43.5, 46, 48.5, 51 cm



38 - 46 cm
39 - 51 cm**




20 - 23 - 26 cm,
24 cm, 19 - 29.5 cm




35 - 48 cm




0° - 6°



7°
0° - 30°




Seat width + 21 cm




90 - 103 cm



105 - 111 cm



83 - 90 cm



30.5 cm



15.2 kg



125 / 160 kg



Transport:
Weight 8,3 kg



Crash tested
ISO 7176-19***

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

ACTION FAMILY CROSS COMPATIBLE ACCESSORIES

	Action ¹ _{NG}	Action ² _{NG}	Action ³ _{NG}	Action ⁴ _{NG}	Action ³ _{Junior}
<p>Fixed 7° backrest, height adjustable push handles (15 cm) :</p> <ul style="list-style-type: none"> • 430 mm height on Adult version • 400 mm height on Junior version 			0425	0314	0425
<p>Reclining backrest from 0° to 30° (ratchet) :</p> <ul style="list-style-type: none"> • 390, 410 and 430 mm heights, delivered in 430 mm on Adult version • 350, 370, 390, 410 and 430 mm heights, delivered in 370 mm on Junior 		●	0435	0320	0390
<p>Reclining backrest from 0° to 30° (gas spring) :</p> <ul style="list-style-type: none"> • 390, 410 and 430 mm, heights, delivered in 430 mm, including backbrace on Adult version • 400 mm height on Junior version 			0445	0333	0395

ACTION FAMILY CROSS COMPATIBLE ACCESSORIES

	Action ¹ _{NG}	Action ² _{NG}	Action ³ _{NG}	Action ⁴ _{NG}	Action ³ _{Junior}
<p>Elevating legrest Junior from 0° to 80°, swing in and out hanger, including calfpad</p> 					0520 0530 0570 0580
<p>Flip up angle and depth adjustable footplate</p> 	0610 0630	●	0610 0630	0490 0510	0605 0625
<p>Amputee legrest</p> 	0545 0595	●	0545 0595	0420 0470	
<p>One piece footplate, angle and depth adjustable</p> 			0611	0491	

ACTION FAMILY CROSS COMPATIBLE ACCESSORIES

		Action ¹ _{NG}	Action ² _{NG}	Action ³ _{NG}	Action ⁴ _{NG}	Action ³ _{Junior}
<p>Neckrest grey Dartex, including backbrace (except for Laguna backrest), adjustable in depth and height</p> <p>NEW</p>		0905	●	0905	0716	
<p>Headrest grey Dartex, including backbrace (except for Laguna backrest), adjustable in depth, width and height</p> <p>NEW</p>				0906	0717	
<p>Removable acrylic tray</p>		1480	●	1480	1270	1470
<p>Removable full length, adjustable to 200, 230, 260 and 280 mm height (n^o1)</p>			●	0650 0750		

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za