

# SATURN PORTERS



<b>NAME</b>	SATURN PORTER WHEELCHAIR
<b>PRODUCT CODE</b>	POR020 / POR020T
<b>SIZE</b>	20"
<b>COLOURS</b>	BLACK
<b>FRAME</b>	NON FOLDING MILD STEEL POWDER COATED REINFORCING PUSH BAR AND BAR BETWEEN FOOTREST HANGERS
<b>ARMS</b>	FIXED FULL ARM
<b>FOOTREST</b>	FIXED FOOTREST WITH REINFORCING BAR
<b>CASTORS</b>	8X2" PU CASTOR ON 2MM STEEL FORK
<b>RIMS</b>	24X1 <sup>3/8</sup> MAG RIM WITH 24X1 <sup>3/8</sup> PU TYRE - FIXED AXLE
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON WITH "Z" FOLD POCKET FOR FILES
<b>STRENGTH</b>	CHROME PUSHRIMS, STEEL SIDE PANELS, RIGID REINFORCED BOX SHAPE DESIGN
<b>WEIGHT</b>	LOAD CAPACITY - 125KG WEIGHT 24KG
<b>PACKING</b>	110 X 80 X 97 (CM)
<b>ACCESSORIES</b>	STAINLESS STEEL DRIP STAND "T" ROD
<b>RIMS</b>	12 1/2 x 2 1/4 MAG W/Solid tyre

# BASIK TRANSIT



<b>NAME</b>	BASIK TRANSIT WHEELCHAIR
<b>PRODUCT CODE</b>	YK9092/3
<b>SIZE</b>	17" & 19"
<b>COLOURS</b>	BLACK
<b>FRAME</b>	MILD STEEL & ALUMINIUM POWDER COATED
<b>ARMS</b>	FIXED FULL ARMS
<b>FOOTREST</b>	REMOVABLE
<b>CASTORS</b>	FOUR 8X1" SOLID CASTOR
<b>RIMS</b>	
<b>UPHOLSTERY</b>	BLACK NYLON
<b>WEIGHT</b>	LOAD CAPACITY – 90KG GROSS WEIGHT – 11KG
<b>PACKING</b>	76 X 26 X 60 (CM)
<b>USE</b>	LIGHT INDOOR / SMOOTH SURFACE USE

# INVACARE ALU-LITE



**Folding Backrest**  
Reduces overall size when folded for transportation or storage.

**Swing away footrest**  
Footplates fold up as shown and footrests swing away, allowing easier access in or out of the chair.

**Foldable cross frame**  
Lightweight frame folds flat for easier portability.

PRODUCT CODE: 1517738

## Technical data

Invacare Alu Lite	405 / 455 mm	430 mm	0°	490 mm	0°	450 mm	80°	340 - 460 mm

Invacare Alu Lite	230 mm	Seat width + 110 mm	250 mm	1030 mm	760 mm	945 mm

Invacare Alu Lite	12,5°	890 mm	100 kg	12 kg	10,2 kg

# INVACARE HEAVY DUTY TRANSPORT CHAIR



<b>NAME</b>	INVACARE HEAVY DUTY TRANSPORT CHAIR
<b>PRODUCT CODE</b>	TRHD22FR
<b>SEAT WIDTH</b>	22"
<b>SEAT DEPTH</b>	18"
<b>BACK HEIGHT</b>	20"
<b>SEAT-TO-FLOOR HEIGHT</b>	19.5"
<b>OVERALL WIDTH OPEN</b>	28"
<b>OVERALL DEPTH W/O RIGGINGS</b>	36.5"
<b>LOAD CAPACITY</b>	400LB

## Features

- 13 month limited warranty on side frames and cross braces
- 12" rear wheels with flat free tyres
- Removable padded desk arms
- Carbon steel frame with dual cross braces
- Swingaway footrest with aluminium footplates
- Flame retardant black nylon upholstery

# YK9070



## YK9070

- Aluminum chair frame
- Flip back Desk Arm
- Tool-free adjustable Swing Away footrest
- 24" spoke wheel with pneumatic tire
- 8" front castor
- Drum brake
- Foldable back for easy storage

ITEM NO.	SIZE	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to seat Height	Weight Capacity (Kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (Kg)
YK9070	18"	103	68	92	33	46	40	40	40	115	8"	24"	16.5/15

# OMEGA FULL RECLINER



<b>NAME</b>	OMEGA FULL RECLINER
<b>PRODUCT CODE</b>	YK9132
<b>SIZE</b>	16", 18", 20"
<b>FRAME</b>	CHROME STEEL
<b>FOOTREST</b>	ELEVATING LEG REST
<b>CASTOR</b>	8X1" SOLID CASTOR
<b>REAR RIM</b>	24X1" MAG WITH 24X1" SOLID TYRE
<b>HEAD REST</b>	REMOVABLE HEAD REST EXTENSION
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON
<b>CARRY CAPACITY</b>	LOAD CAPACITY – 100KG
	WEIGHT – 22KG
<b>ARMS</b>	REMOVABLE FULL ARM WITH UPHOLSTERED ARM PAD
<b>PACKING</b>	31X92X96CM

# INVACARE 9000XT FULL RECLINER



Invacare® 9000 XT Recliner  
Wheelchair Features  
Dynamic recline range from 95°  
to 175°, +/- 5 °  
Rear axle base moves,  
providing stable centre-of-gravity  
in recline position.  
Multiple wheel offerings allow  
variation of seat-to-floor heights  
Urethane read tires mounted on  
“no flex” wheels offer superior  
performance  
Urethane casters for better ride  
and lighter weight  
Low-maintenance, dual-sealed  
precision bearings  
Heavy-duty inner liners keep  
seat and back from stretching  
Stocked with attendant  
(standard) or patient-operated  
wheel locks

<b>NAME</b>	INVACARE 9000XT FULL RECLINER
<b>PRODUCT CODE</b>	9000XT
<b>SIZE</b>	14, 16, 18, 20, 22"
<b>FRAME</b>	MILD STEEL BLACK PRISM POWDER COATED
<b>FOOTREST</b>	REMOVABLE
<b>CASTOR</b>	6" / 8" URETHANE CASTORS
<b>REAR RIM</b>	URETHANE REAR TIRES
<b>HEAD REST</b>	10" REMOVABLE HEADREST
<b>UPHOLSTERY</b>	BLACK NYLON
<b>CARRY CAPACITY</b>	LOAD CAPACITY 250lb. WEIGHT - 39lb/18kg
<b>ARMS</b>	FLIP BACK FIXED LENGTH
<b>PACKING</b>	32 X 13 X 37"

# INVACARE ACTION 3 COMFORT



<b>NAME</b>	INVACARE ACTION 3 COMFORT
<b>PRODUCT CODE</b>	CEG0041
<b>SIZE</b>	16", 18", 20"
<b>FRAME</b>	FOLDING POWDER COATED ALUMINIUM COLOUR: RED
<b>BACKREST</b>	RECLINING WITH GAS ASSISTED SHOCKS
<b>FOOTREST</b>	REMOVABLE ELEVATING
<b>CASTOR</b>	8X1" SOLID CASTOR
<b>REAR RIM</b>	24X13/8" ALUMINIUM SPOKE RIM WITH 24X13/8" GREY PU TYRE/REMOVABLE
<b>UPHOLSTERY</b>	GREY DARTEX COMFORT SEAT & BACK WITH HEAD REST
<b>CARRY CAPACITY</b>	LOAD CAPACITY – 125KG WEIGHT – 21KG
<b>ARMS</b>	HEIGHT ADJUSTABLE COMFORT ARMS
<b>PACKING</b>	108X92X42CM



# INVACARE CLEMATIS



<b>NAME</b>	INVACARE REA CLEMATIS
<b>PRODUCT CODE</b>	156087
<b>SIZE</b>	96 -112 CM
<b>FRAME</b>	PEARL GREY
<b>BACKREST</b>	REMOVABLE BACKREST
<b>FOOTREST</b>	ADJUSTABLE BODY LENGTH LEGRESTS
<b>CASTOR</b>	8 X 2" SOLID
<b>REAR RIM</b>	24 X 1 <sup>3/8</sup> ALUMINIUM RIM REMOVABLE
<b>UPHOLSTERY</b>	GREY DARTEX WITH REMOVABLE HEAD REST
<b>CARRY CAPACITY</b>	LOAD CAPACITY 125KG
<b>ARMS</b>	WEIGHT - 20,5 - 21,5 KG ARMREST WITH ANGLE
<b>PACKING</b>	H 59,5 x L 84 x 68,5 cm

The Rea™ Clematis® complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information vehicular transportation for this wheelchair, please consult the following pages on our website: [www.invacare-rea.com](http://www.invacare-rea.com) under Rea wheelchairs/Technical downloads.

Width 39 (h 59.5 L 84 x 58.5 cm)  
Width 44 (h 59.5 x L 84 x 63.5 cm)  
Width 49 (h 59.5 x L 84 x 68.5 cm)

# MATRIX ELITE BACK



PRODUCT CODE: PBE

**Invacare® Matrix® Elite Back**  
Model no PBE

### Features and Function Benefits

#### Outer Cover

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available



3" CONTOUR DEPTH



SACRAL /  
LUMBAR  
SUPPORT

# MATRIX ELITE DEEP BACK



PRODUCT CODE: EDB

**Invacare® Matrix® Elite Deep Back** Model No: EDB

### Features and Functional Benefits

#### Outer cover:

- Oversized cover/foam for extra comfort and protection
- Moisture resistant and breathable
- Includes removable lumbar positioning pad
- Optional privacy flap available
- Hardware:
  - "Easy Set" hardware installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair
  - Refer to hardware page
- Back shell:
  - Stylish and breathable back shell
  - Lightweight – most models under 4lb
- Sacral / Lumbar support:
  - Includes removable positioning pad

# MATRIX HARDWARE



**Hardware:**

- "Easy Set" hardware installed and adjusted in minutes – adjustments can be made while user is seated in the wheelchair
- Refer to hardware pageBack shell
- Stylish and breathable back shell
- Lightweight – most models under 3lb
- Sacral/Lumbar support
- Includes removable positioning pad

► SWING-AWAY	► OFFSET SWING-AWAY	► LATERAL PADS	► FIXED	► OFFSET FIXED
		<p>Extra Small</p> <p>Small</p> <p>Medium</p> <p>Large</p>		
<p>► 2 1/2" of vertical adjustment on back pan</p> <p>► Mount inside back shell without interfering w/depth/angle adjustment</p> <p>► Easy push button release</p> <p>► Not offered on Elite Deep</p> <p>► 4 pad sizes: XS, S, M &amp; L</p>	<p>► Provides an additional 3" of vertical position (higher or lower) relative to mounting slots</p> <p>► Mount inside back shell without interfering w/depth/angle adjustment</p> <p>► Easy push button release</p> <p>► Not offered on Elite Deep</p> <p>► 4 pad sizes: XS, S, M &amp; L</p>	<p>► Lateral Pads are available in 4 pad sizes: XS, S, M &amp; L</p> <p>XS: 3.25" high x 4.25" long</p> <p>S: 3.5" high x 5.25" long</p> <p>M: 3.75" high x 5.75" long</p> <p>L: 4.25" high x 6.5" long</p>	<p>► 2 1/2" of vertical adjustment on back pan</p> <p>► Mount inside back shell without interfering with depth/angle adjustment</p> <p>► Available in standard and pediatric sizes</p> <p>► Not offered on Elite Deep</p>	<p>► Provides an additional 3" of vertical position (higher or lower) relative to mounting slots</p> <p>► Not offered on Elite Deep</p>

# Medop<sup>cc</sup>



Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

