

TIAN LIGHTWEIGHT FOLDING POWER WHEELCHAIR

CODE: DYW-459-46A8



- Aluminium powder coating frame, can be folded in small size;
- Flip-up armrest; foldable footplate;
- With pressure-reducing, washable, detachable “S” shape upholstery that are made by breathable material;
- Forward and backward adjustable controller can be installed in left armrest or right armrest freely;
- Motor: Brushless 180W x 2 with Electromagnetic Brake;
- Range: 20 ± 10% km;
- Max. speed: 6km/h;
- 8*PU casters; 8” sold rear wheel (Drive wheel);
- N.W.:26.4KGs lithium battery included;
- Max. user weight: 100KGs;
- Lithium battery: 24V/10Ah;

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

BASIK STANDING

PRODUCT CODE ALK158



Model /**ALK158**
Name /*Electric Wheelchair*



- 1) Spray powder frame with breathable mesh cushion
- 2) 8" big castor, 12" PU rear wheel
- 3) Detachable plate, the armrest can be turned backwards
- 4) can go up and down and walk upright by electric motor
- 5) with safe light
- 6) Seat width of 40cm

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za