

# ALBER E-MOTION



## Power assist for wheelchairs

A manual wheelchair makes it possible to move around and promotes activity. For most wheelchair users it is often impossible to access all areas of everyday life - even with adequate physical strength. Slopes, kerbstones, thresholds regularly present insurmountable obstacles. The new **Alber e-motion** power assist allows wheelchair drivers to keep and even raise their level of active mobility, and provides a variety of therapeutic benefits.

### Technical data

<b>e-motion</b>	 max. 6 km/h	 25 km (according to ISO 7176)	 2 x 60 W	 1 kg	 22 kg (one wheel: 11 kg)	 130 kg
-----------------	--	--	---	---	--	---

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za