

KUSCHALL CHAMPION SK

Sleek, smart styling meets compactness and practicality. The active rigid wheelchair that folds is now even smaller.

Having the benefit of an everyday wheelchair with such small folded dimensions makes light work of all aspects of travel. It's the perfect companion in your busy day to day commuting and essential for any holiday destination.

Evolved from our highly popular Champion, the SK features a folding front frame, the result is a wheelchair that is lightweight, easy to handle and convenient to transport.


Exceptional driving performance, rigidity and quality are also key features, ensuring as well as being practical you can be sure it delivers the strength and durability needed for an active lifestyle.





THE MOST FOLDABLE - SMALL BUT SMART


TECHNICAL DATA


 **Frame angle**
75°/85°
Aluminium folding front frame

 **Seat width**
360 – 460 mm
in increments of 20 mm

 **Seat depth**
340 – 480 mm
in increments of 20 mm

 **Seat height**
front: 450 – 540 mm
rear: 390 – 490 mm
stepless adjustable


 **Knee-to-heel length**
380 – 500 mm
stepless adjustable

 **Backrest height**
300 – 465 mm
in increments of 15 mm


 **Backrest angle**
76°/80.5°/85°/89.5°/94°


 **Max. safe slope**
7°


 **Reversing width**
920 – 1230 mm


 **Total width**
SW + 160 mm

 **Total folded dimensions**
Rear 280mm
Front 330mm
Length 670mm
Width 220mm

 **Total length**
75°: approx. 930 mm
85°: approx. 850 mm

 **Total weight unloaded**
approx. 9.2 kg

 **Transport weight w/out rear wheels**
approx. 6.7 kg

 **Max. user weight**
120 kg

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za