

# KUSCHALL K-SERIES ATTRACT





## FEATURES AND BENEFITS


- 1. CONFIGURABLE:** It is available in a wide range of seat sizes and seat heights. Combined with a wide range of options, accessories and colours, the K-Series can be tailored to suit your exact needs and personality.
- 2. ADJUSTABLE:** Adjustments come as standard, so finding your perfect position is guaranteed!  
Possible Adjustments:  
Front seat height                      Rear seat height  
Centre of gravity                      Backrest height and angle
- 3. EASY TRANSPORT:** Quick and easy folding with the perfect open frame design ensures effortless car transfers. Keep independent and on the go!
- 4. MATRX BACKRESTS AND CUSHIONS:** Our range offers additional posture positioning giving you the freedom to enjoy the benefits of life in an active wheelchair.

AN ATTRACTIVE WHEELCHAIR FOR EVERYONE.


## TECHNICAL DATA

 **Frame angle**  
75°/90°  
V-frame or tapered


 **Seat width**  
320 – 500 mm  
in increments of 20 mm


 **Seat depth**  
350 – 525 mm  
in increments of 25 mm

 **Seat height**  
front: 450 – 520 mm  
rear: 380 – 490 mm  
in increments of 10 mm


 **Knee-to-heel length**  
290 – 480 mm  
in increments of 10 mm


 **Backrest height**  
270 – 480 mm  
in increments of 15 mm


 **Backrest angle**  
74°/78°/82°/86°/90°


 **Max. safe slope**  
7°

 **Reversing width**  
920 – 1230 mm


 **Total width**  
SW + 170 mm

 **Centre of gravity**  
80 mm/95 mm/110mm/  
125 mm/140 mm

 **Total length**  
75°: approx. 830 - 1055 mm  
90°: approx. 760 - 985 mm

 **Total weight unloaded**  
approx. 7.8 kg

 **Transport weight w/out rear wheels**  
approx. 5.2 kg

 **Max. user weight**  
130 kg

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za