

KUSCHALL K-SERIES

PRODUCT CODE: KUSK

It was, everything with the **K-Series** was incredibly easy, unique and light. With its broad selection of possible configurations, different materials and the wide range of colours to choose from, the **K-Series** brings a new dimension to everything you need.

And not only is it more versatile. It is a wheelchair that gets you on the move much more easily and efficiently and on an entirely new scale.



Kuschall wheelchairs are compatible with Matrix backs and cushions



A «no limits» range



Lightness and performance



Backrest adapter plate with lasered swooshes®



Front seat height adjustable



TECHNICAL DATA K-SERIES



Frame angle
75° / 90°
V-frame or tapered



Knee-to-heel length
290 – 480 mm
in increments of 10 mm



Total width
SW + 170 mm



Total length
75°: approx. 830 mm
90°: approx. 760 mm



Seat height
front: 450 – 520 mm
rear: 380 – 490 mm
in increments of 10 mm



Reversing width
920 – 1230 mm



Max. user weight
130 kg



Crash Test approved
Fulfills the requirements
of ISO 7176-19



Max. safe slope
7°



**Transport weight
w/out rear wheels**
approx. 5.1 kg



Backrest height
270 – 480 mm
in increments of 15 mm



**Total weight
unloaded**
approx. 7.8 kg



Seat depth
350 – 525 mm
in increments of 20 mm



Backrest angle
74° / 78° / 82° / 86° / 90°



Seat width
320 – 500 mm
in increments of 20 mm



Centre of gravity
80 mm / 95 mm / 110mm /
125 mm / 140 mm

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za