

# KUSCHALL ULTRA LIGHT

PRODUCT CODE: KUSUT













## Options



These are only a few of the many options available for this model. Contact your nearest Kuschall specialist for more information.



## Technical Data


 Seat width: 28 – 50 cm (in increments of 2 cm)	 Backrest angle: 82°/86°/90°/94°/98°/102°
 Seat depth: 32 – 50 cm (in increments of 2 cm)	 Total width (with 3° camber): SW + 17 cm total width folded approx. 29 cm
 Seat height: front: 43 – 53 cm rear: 39 – 50 cm (in increments of 1 cm)	 Total length: 70°: 88 cm 80°: 83 cm
 Backrest height: 30 – 51 cm (in increments of 1.5 cm)	 Total weight unloaded: Approx. 9,9 kg
 Knee-to-heel length: 20 – 50 cm (in increments of 1 cm)	 Max. user weight: 130 kg
	 Transport weight w/out rear wheels, Approx. 7,5 kg

All measurements and weights indicated are based on seat width 40 cm and seat depth 40 cm in the lightest and most compact configuration possible. Weight and measurements can change with the different configurations.

## Colours



## CPS (Colour Power Set)



# Medop cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za