

INVACARE A4



Frog Legs® front suspension forks absorb initial impact slowing vibration before it gets to the user.



1. Undermount wheel locks
2. Fold-down back



LSS dual rear shocks absorb energy independently for added stability.

Adjustability for a Great Fit

It is critical for rigid chairs consumers to have a chair that can be adjusted to match their level of ability and support their active lifestyles. The A-4 Wheelchairs are a great fit for these demanding consumers. These chairs can be adjusted for maximum stability or fine tuned for peak performance – quickly and easily, today and in the future.

Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Back Angle	Front Frame Angle	Camber	Tubing Diameter
TA4	Seat Width + 6.75" - 14.25"	12"-22" W 14"-22" D	17"-21" Front 14.5"- 20" Rear	85°, 90°, 95°, 100°, 105°	80°	0°, 3°, 6°, 9°, 12°	1"
	Knee-to-Heel Range	Center-of-Gravity	Wheels	Casters	Total Weight*	Weight Capacity	Limited Warranty
	11.25"-20.25" 1" Increments	Adjustable	22", 24", 25", 26"	3", 4", 5", 6"	14.5 lb.	300 lb.	5 Years Frame

* 15"x16" without rear wheels, folding back, adjustable tension upholstery and 4" casters.

Please see the Owner's Manual or Invacare Configurator for complete option/specification combinations.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za