

GRIFFIN SEMI-ELECTRIC STAND-UP



| | |
|--------------------------|--|
| NAME | SEMI-ELECTRIC STAND-UP |
| PRODUCT CODE | 50003 |
| CONTROLS | SINGLE ACTUATOR, FINGER OPERATED |
| CASTORS | 6X1" GREY CASTORS |
| REAR RIM | 24X1 3/8" PU GREY TYRE ON ALUMINIUM SPOKE RIM |
| BATTERIES | 2 X SEALED 3.8AH/24VOLT |
| LOAD CAPACITY | LOAD CAPACITY - 110KG WEIGHT - 30 KG |
| ADJUSTABILITY | BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY |
| PACKING SIZES | 65X68X94CM 14", 16", 18" |

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za