

# INVACARE ACTION 2

## The Next Generation has arrived!


The all-new **Action<sup>2</sup><sub>NG</sub>** has been developed from the tried-and-tested **Action<sup>3</sup>** concept. The specifications meet the basic needs of the users whilst refusing to compromise on safety or comfort. The next generation **Action<sup>2</sup><sub>NG</sub>** shares the same essential platform as the **Action<sup>3</sup>**, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the **Action<sup>2</sup><sub>NG</sub>** perfectly fits into the Action family.



### Features



**Swing in and out footrest**  
The 80 degree footrest swings in and out for easy transfers in or out of the chair.



**Adjustable center of gravity**  
The newly-designed rear wheel plate can be adjusted by 5 cm.



**Anti-tippers**  
Provide extra safety on slopes or uneven ground.

### Technical data

38, 40,5, 43, 45,5, 48 cm	40,5 / 43 cm	48,5, 51 cm	43 cm	20 - 26 cm	35 - 48 cm	80° / 0 - 80°	7° / 0° - 30°
SW+ 20,5 cm	93,5 cm	82 / 84,5 cm	41 cm	15 kg	max. 125 kg	Crash tested ISO 7176-19*	

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za