

INVACARE ACTION 4

Where stability meets style!

A comfortable, foldable, lightweight wheelchair, the *Invacare® Action⁴_{NG}* is ideally suited for everyday use.

A key feature of this next-generation model is the innovative Invacare® hallmark of cross-compatibility; a shared platform with other models in the Action family as well as other wheelchairs in our range means that users and advisors can easily update or refurbish their wheelchairs and enjoy maximum opportunities for modularity and personalized adjustment.



Options

NEW



Heavy Duty Version 160 kg
A wider choice of seat widths (55.5 & 60.5 cm), reinforced chassis and a new greater carrying capacity (maximum user weight 160 kg).




Hemi version*
A 38.5 cm seat height, dual handrim propulsion, hemi armrest and built out castor housing.



Amputee Legrest
The amputee legrest is both height- and angle- adjustable.

Technical data




38, 40.5, 43, 45.5, 48, 50.5, 55.5, 60.5 cm




40, 45 and 50 cm



38.5*, 41, 43.5, 46, 48.5, 51 cm



38 - 46 cm
39 - 51 cm**




20 - 23 - 26 cm,
24 cm, 19 - 29.5 cm




35 - 48 cm



0° - 6°



7°
0° - 30°



Seat width + 21 cm




90 - 103 cm




105 - 111 cm



83 - 90 cm



30.5 cm



15.2 kg



125 / 160 kg



Transport:
Weight 8,3 kg



Crash tested
ISO 7176-19***

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za