




ACTION FAMILY CROSS COMPATIBLE ACCESSORIES

	Action ¹ _{NG}	Action ² _{NG}	Action ³ _{NG}	Action ⁴ _{NG}	Action ³ _{Junior}
<p>Fixed 7° backrest, height adjustable push handles (15 cm) :</p> <ul style="list-style-type: none"> • 430 mm height on Adult version • 400 mm height on Junior version 			0425	0314	0425
<p>Reclining backrest from 0° to 30° (ratchet) :</p> <ul style="list-style-type: none"> • 390, 410 and 430 mm heights, delivered in 430 mm on Adult version • 350, 370, 390, 410 and 430 mm heights, delivered in 370 mm on Junior 		●	0435	0320	0390
<p>Reclining backrest from 0° to 30° (gas spring) :</p> <ul style="list-style-type: none"> • 390, 410 and 430 mm heights, delivered in 430 mm, including backbrace on Adult version • 400 mm height on Junior version 			0445	0333	0395