

# INVACARE Action Family: The changing face of Mobility

**The Action Family product offer consists of:**

Action2 - Standard aluminium wheelchair  
Action3 Junior - Paediatric wheelchair from age 3 and up

Action3 - Lightweight wheelchair, single crossbar, highly configurable with a huge range of options, including:

- Hemiplegic version
- One arm drive version
- Double hand rim version
- Comfort version
- and more

Action4 - Lightweight wheelchair, double crossbar, highly configurable with huge range of options, including:

- Hemiplegic version
- One arm drive version
- Double hand rim version
- Comfort version

Heavy-duty option up to 160kg user weight



Cross Compatibility table: Action<sup>2</sup> Action<sup>3</sup>

Leg rest	x	x
Footplates	x	x
Armrest	x	x
Arm pads	x	x
Side guards	x	x
Seat upholstery	x	x
Backrest canes	x	x
Backrest upholstery		x
Push handles	x	x
Brakes	x	x
Fork		x
Castor wheels	x	x
Frame	x	x
User weight	125KG	125KG

Action<sup>3</sup> One-arm drive version



Action<sup>3</sup> Comfort version



Action<sup>3</sup> Hemiplegic version



# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za