

YK9070



YK9070

- Aluminum chair frame
- Flip back Desk Arm
- Tool-free adjustable Swing Away footrest
- 24" spoke wheel with pneumatic tire
- 8" front castor
- Drum brake
- Foldable back for easy storage

ITEM NO.	SIZE	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to seat Height	Weight Capacity (Kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (Kg)
YK9070	18"	103	68	92	33	46	40	40	40	115	8"	24"	16.5/15