

# INVACARE 9000XT FULL RECLINER



## Invacare® 9000 XT Recliner Wheelchair Features

Dynamic recline range from 95° to 175°, +/- 5°  
Rear axle base moves, providing stable centre-of-gravity in recline position.  
Multiple wheel offerings allow variation of seat-to-floor heights  
Urethane tread tires mounted on “no flex” wheels offer superior performance  
Urethane casters for better ride and lighter weight  
Low-maintenance, dual-sealed precision bearings  
Heavy-duty inner liners keep seat and back from stretching  
Stocked with attendant (standard) or patient-operated wheel locks

<b>NAME</b>	INVACARE 900XT FULL RECLINER
<b>PRODUCT CODE</b>	9000XT
<b>SIZE</b>	14, 16, 18, 20, 22"
<b>FRAME</b>	MILD STEEL BLACK PRISM POWDER COATED
<b>FOOTREST</b>	REMOVABLE
<b>CASTOR</b>	6" / 8" URETHANE CASTORS
<b>REAR RIM</b>	URETHANE REAR TIRES
<b>HEAD REST</b>	10" REMOVABLE HEADREST
<b>UPHOLSTERY</b>	BLACK NYLON
<b>CARRY CAPACITY</b>	LOAD CAPACITY 250lb. WEIGHT - 39lb/18kg
<b>ARMS</b>	FLIP BACK FIXED LENGTH
<b>PACKING</b>	32 X 13 X 37"

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za