

# OMEGA FULL RECLINER



<b>NAME</b>	OMEGA FULL RECLINER
<b>PRODUCT CODE</b>	YK9132
<b>SIZE</b>	16", 18", 20"
<b>FRAME</b>	CHROME STEEL
<b>FOOTREST</b>	ELEVATING LEG REST
<b>CASTOR</b>	8X1" SOLID CASTOR
<b>REAR RIM</b>	24X1" MAG WITH 24X1" SOLID TYRE
<b>HEAD REST</b>	REMOVABLE HEAD REST EXTENSION
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON
<b>CARRY CAPACITY</b>	LOAD CAPACITY – 100KG WEIGHT – 22KG
<b>ARMS</b>	REMOVABLE FULL ARM WITH UPHOLSTERED ARM PAD
<b>PACKING</b>	31X92X96CM

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za