

# SATURN PORTERS



removable "T" bar



Available with 12" rear rims option

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| <b>NAME</b>         | SATURN PORTER WHEELCHAIR                                                                         |
| <b>PRODUCT CODE</b> | POR020 / POR020T                                                                                 |
| <b>SIZE</b>         | 20"                                                                                              |
| <b>COLOURS</b>      | BLACK                                                                                            |
| <b>FRAME</b>        | NON FOLDING MILD STEEL POWDER COATED<br>REINFORCING PUSH BAR AND BAR BETWEEN<br>FOOTREST HANGERS |
| <b>ARMS</b>         | FIXED FULL ARM                                                                                   |
| <b>FOOTREST</b>     | FIXED FOOTREST WITH REINFORCING BAR                                                              |
| <b>CASTORS</b>      | 8X2" PU CASTOR ON 2MM STEEL FORK                                                                 |
| <b>RIMS</b>         | 24X1 <sup>3/8</sup> MAG RIM WITH 24X1 <sup>3/8</sup> PU TYRE - FIXED<br>AXLE                     |
| <b>UPHOLSTERY</b>   | BLACK REINFORCED NYLON WITH "Z" FOLD<br>POCKET FOR FILES                                         |
| <b>STRENGTH</b>     | CHROME PUSHRIMS, STEEL SIDE PANELS, RIGID<br>REINFORCED BOX SHAPE DESIGN                         |
| <b>WEIGHT</b>       | LOAD CAPACITY - 125KG WEIGHT 24KG                                                                |
| <b>PACKING</b>      | 110 X 80 X97 (CM)                                                                                |
| <b>ACCESSORIES</b>  | STAINLESS STEEL DRIP STAND "T" ROD                                                               |
| <b>RIMS</b>         | 12 ½ x 2 ¼ MAG W/Solid tyre                                                                      |

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za