

ACTION 3 JUNIOR EVOLUTIF

A wheelchair that grows with the child

Invacare **Action³ Junior** is a lightweight, foldable paediatric wheelchair designed for children aged between 3 and 15 years. **Action³ Junior** has been developed to match the individual needs of the child and can grow as they grow. Because a child's needs will change the **Action³ Junior** offers a large range of options to accompany the child in their development and to provide the necessary clinical support. The **Action³ Junior** comes in attractive, modern styles, designed to appeal to kids.



Seat-width adjustment
The seat width of the chassis can be adjusted with an innovative telescopic crossbar for four fixed-seat positions. (Child: 20,5 to 28 cm, Young: 30,5 to 38 cm.)



Upholstery adjustment
As the chassis is made wider, so the upholstery on the back and on the seat of the chair can be made wider with Velcro flaps.



Table tray
In addition to the standard slide-on table tray, we offer a swing-away tray which can easily be folded away.



Adjustable seat height and balance
For improved driving performance and different height options, the rear wheel can be mounted in 20 different positions.



Technical data

Invacare Action³ Junior Child	20,5 / 23 / 25,5 / 28 cm	20 - 25 / 25 - 30 / 30 - 35 cm	41 / 43,5 / 46 cm	37 / 40 / 43 cm ¹	15 - 25 ² / 18 - 28 cm ³ 13 - 27 ⁴ / 22 - 32 ⁵ / 30 - 45 cm ⁶	
Invacare Action³ Junior Young	30,5 / 33 / 35,5 / 38 cm	27,5 - 32,5 / 32,5 - 37,5 cm	41 / 43,5 / 46 cm	37 / 40 / 43 cm ¹	15 - 25 ² / 18 - 28 cm ³ 13 - 27 ⁴ / 22 - 32 ⁵ / 30 - 45 cm ⁶	
Invacare Action³ Junior Child	3° / ±1°	7° / 0° - 30° ⁷	SW + 19 / + 17,5 cm ⁸	87 - 113 cm	67 - 97 cm	32,5 cm
Invacare Action³ Junior Young	3° / ±1°	7° / 0° - 30° ⁷	SW + 19 / + 17,5 cm ⁸	87 - 113 cm	67 - 97 cm	35 cm
Invacare Action³ Junior Child	12,3 kg	80 kg	7,5 kg			
Invacare Action³ Junior Young	15 kg	80 kg	9 kg			

1/ Look in the EPF for possible adjustment of upholstery and pushhandles

2/ Desk

3/ Junior

4/ 90 degree hanger

5/ 80 degree short hanger

6/ 80 degree long hanger

7/ Recliner

8/ Transit Version

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za