

INVACARE ACTION 3 JUNIOR



NAME	INVACARE ACTION 3 JUNIOR
PRODUCT CODE	CNB0030
SIZE	CHILD 8"-1" JUNIOR 12"-15"
COLOURS	OPTIONAL
FRAME	FOLDING ALUMINIUM FRAME
ARMS	REMOVABLE & HEIGHT ADJUSTABLE "T" ARM
FOOTREST	REMOVABLE
CASTORS	6" SOLID CASTER
RIMS	20" ALUMINIUM SPOKE RIM WITH PU TYRE
	QUICK RELEASE AXLES
UPHOLSTERY	BLACK REINFORCED NYLON
ADJUSTMENT	CENTRE OF GRAVITY
	LOW SLUNG
WEIGHT	LOAD CAPACITY – 80KG
	GROSS WEIGHT – 15KG

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za