

## SATURN ALU-LITE KIDDIES



<b>NAME</b>	SATURN ALU-LIGHT
<b>PRODUCT CODE</b>	LW25-77K
<b>SIZE</b>	14" X 14", 14" X 16"
<b>COLOURS</b>	METALLIC BLUE, WET BLACK
<b>FRAME</b>	FOLDING ALUMINIUM POWDER COATED
<b>ARMS</b>	FIXED HEIGHT FLIP BACK FULL ARM
<b>FOOTREST</b>	REMOVABLE
<b>CASTORS</b>	8X1" SOLID CASTOR ON ALUMINIUM FORK
<b>RIMS</b>	24X13/8 MAG RIM WITH PU TYRE 12MM QUICK RELEASE AXLES
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON
<b>ADJUSTMENT</b>	CENTRE OF GRAVITY, LOW SLUNG HEIGHT ADJUSTABLE BACK REST WITH PUSH HANDLES
<b>WEIGHT</b>	LOAD CAPACITY – 115KG GROSS WEIGHT – 17KG (100X34X85CM)

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za