

SATURN KIDDIES



NAME	SATURN KIDDIES
PRODUCT CODE	KIDDIES
SIZE	12"
COLOURS	WET BLACK
FRAME	FOLDING MILD STEEL POWDER COATED
ARMS	HEIGHT ADJUSTABLE FULL ARM
FOOTREST	REMOVABLE 90° FOOTREST
CASTORS	8X2" PU CASTOR ON STEEL FORK
RIMS	20 X 1.75 STEEL SPOKE RIM WITH 24X1.75 TYRE WITH NO MORE FLAT (NMF), FIXED AXLES
UPHOLSTERY	BLACK REINFORCED NYLON
WEIGHT	LOAD CAPACITY – 50KG GROSS WEIGHT – 16KG
PACKING	100X34X85CM

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za