


# FREEDOM DESIGN NXT



## Transfers as easy as A, B, C...

- A:** Travels light, weighing just 22 pounds.
- B:** Folds into small package.
- C:** Fits into small spaces.



### Seat Frame Widths

The NXT is available in 11 sizes from 10" wide to 20" wide.  
 Each size from 10" wide to 16" wide adjusts an additional 2" in width.  
 The 17" wide size adjusts an additional 1" in width to 18" wide.  
 The 18", 19" and 20" wide sizes are set.

### Seat Frame Depths

The NXT is available in 3 sizes: a short, a medium or a long.  
**NEW NXT Mini:**  
 For x-short 8" to 12" range.  
 Each size adjusts an additional 4" in depth.  
 Seat frame depth is measured without seating, from the front of the back post to the front end of the seat frame.  
 The short range is 12" to 16."  
 The medium range is 14" to 18."  
 The long range is 16" to 20."

### Seat Depth Extensions

For additional seat frame depths, extensions can be added to the back of the seat frame. There is a 2", a 3" or a 4" extension to choose from.

### Seat to Floor Height

The NXT can be configured in a variety of seat to floor heights.  
 Selections in front casters and rear wheel sizes determine a range of choices from 12.5" low to 18" high, and most everything in between.

### Rider Capacity

225 lbs.

### Frame

Available in 25 colors.

### Tilt Range

45 degrees.

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za