

DRAGON JUNIOR POWER WHEELCHAIR

The *Invacare Dragon Junior* is a powerful and adjustable power wheelchair to meet everyday needs of the child. The multiple seat adjustments and the motor options make this wheelchair easy to use.



Technical datas

Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Seat angle	Backrest angle
Dragon Junior 290-360 mm 320-400 mm	300 - 360 mm	425 - 475 mm	450 mm	180 - 250 mm	200 - 350 mm	0° à 17° (powered or man.)	From -10° to 30° (man.) or from 5 to 30° (powered)
Driving unit width	Total height	Total length inc. footrests	Total length without footrests	Total weight	Max. user weight	Battery capacity	Turning diameter
Dragon Junior 610 mm	950 mm	1030 mm	850 mm	100 kg	90 kg	2 x 60 Ah	1660 mm
Max. Climbable obstacle height	Max. safe slope	Motor capacity	Driving range	Speed			
Dragon Junior 60 mm	18 %	2 x 180 W	30 km	6 km/h			

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za