

ALPHA A6 POWER WITH TILT



| | |
|------------------------------|--|
| NAME | ALPHA TILT POWER |
| PRODUCT CODE | A6POW |
| SIZE | 14", 16", 18" |
| BATTERIES | 2 X 36AH/12 VOLT |
| MOTORS | 250 WATT |
| CASTORS | 8 X 2" FOAM FILLED GREY |
| SPEED | 6KM/H |
| RANGE | 20KM (TERRAIN AND USER DEPENDENT) |
| REAR RIMS & TYRES | 12 1/2 X 2 1/4 ALUMINIUM RIM 12 1/2 X 2 1/4 ALUMINIUM RIM |
| ARMS | HEIGHT ADJUSTABLE T ARM |
| GRADIENT | 12 DEGREES |
| CONTROLS | DYNAMIC SHARK REMD11 POWER MODULE PMB11 |
| SEAT FUNCTION | |
| LOAD CAPACITY | LOAD CAPACITY – 120KG WEIGHT – 70KG |
| PACKING | 70X70X100CM |

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za