

DELTA A7 POWER WITH SEAT ELEVATE



NAME	DELTA A7 POWER WITH SEAT RAISER
PRODUCT CODE	A7POW
SIZE	18"
BATTERIES	2X36AH/12 VOLT
MOTORS	25 WATT
CASTORS	8X2" FOAM FILLED GREY
SPEED	6KM/H
RANGE	20KM (TERRAIN AND USER DEPENDENT)
REAR RIMS & TYRES	121/2 X 21/4 ALUMINIUM RIM 121/2 X 21/4 ALUMINIUM RIM
ARMS	HEIGHT & WIDTH ADJUSTABLE ARM
GRADIENT	12 DEGREES
CONTROLS	DYNAMIC SHARK REMD11 POWER MODULE PMB11
SEAT FUNCTION	SEAT ELEVATE MODEL
LOAD CAPACITY	LOAD CAPACITY - 115KG WEIGHT - 85KG
PACKING	71X74X101CM

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za