

INVACARE DRAGON VERTIC

Access all areas – with the push on a button

The Invacare® Dragon Vertic is an agile power chair that not only provides mobility to an active user, but will bring you face to face with other people. Its comfortable seating system easily moves you into a standing position by a simple joystick command, providing increased access and excellent mobility.



NAME	INVACARE DRAGON VERTIC
PRODUCT CODE	VERTIC
CONTROLS	SHARK II STEERING
CASTORS	MINI CATORS
REAR RIM	PUNCTURE PROOF REAR TYRES
STD ACCESSORIES	
BATTERIES	2X12V 60aH
LOAD CAPACITY	LOAD CAPACITY – 131KG WEIGHT – 100KG
ADJUSTABILITY	ARM REST, HEAD REST, SLIDING BACK

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za