

# MERITS P182 JUPITER POWER WHEELCHAIR

## Heavy Duty Folding Power Wheelchair



### Heavy Duty

- Reinforced frame with double cross braces provide a sturdy support for heavy users
- Auxiliary hand brakes offer extra safety
- Four poles motor provides more power and torque
- Optional pneumatic casters and drive wheels

### Adjustable features

- Patented folding battery tray allows improved portability
- Programmable controller with various choices of different performance to fit various needs of users
- Height adjustable flip-up desk armrests and swing away footrests for easier access
- Adjustable mounting bracket for the controller, swing back mounting bracket optional




Four poles motor



### WEIGHT (kg/lb)

capacity	tvl. weight
272/600	83/183

### DIMENSIONS (cm/in)



### TECHNICAL SPECS

Weight capacity: 272kg / 600lbs  
 Max. speed: 6.4kph / 4mph  
 Range: 40km / 25 mi.  
 Turn radius: 66cm / 26"  
 Ground clearance: 12cm / 4.8"

Motor: DC 24V, 300W, 4P(Rating) / 1450W(Max)  
 Controller: PG VSI 70  
 Battery: 12V/ 22NF x 2pcs  
 Charger: 5A off-board  
 Gradient: 8°(600lbs), 10°(525lbs)

Caster: 9"(2.80/2.50-4) foam filled tire  
 Drive wheel: 12-1/2" X 2-1/4" foam filled tire  
 Brake: Intelligent, regenerative, electromagnetic brakes

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za