

MERITS P201 NEPTUNE POWER

Multi- Function Rehab Power Wheelchair

Smart design

- Various seating positions, tilt-in-space seat and recline back, both with manual and power drive
- Modular construction for quick dismount and simple maintenance
- Swing in and out footrests provide easy access to objects

Adjustable fit

- Armrest height adjusted range from 8" to 12", seat width adjustment is 2" from each side
- Electric seat tilt adjustment from 3.5°-16.5°
- Electric back recline adjustment from -5°-45°
- Detachable armrest and legrest for ease of transfer
- Optional power ELR, angle adjustable from 10°-85°



Seat Tilt: Electric seat tilt adjustment from 3.5°-16.5°



Back Recline: Electric back recline adjustment from 85°-135°



Full Recline



WEIGHT (kg/lb)	ttl. w/batt.			
capacity	16"	17"	18"	20"
136/300	91/200	92/202	93/204	95/209

DIMENSIONS (cm/in)

	40/16	40/16	48/19	51-53/20-21	41-53/16-21	16" 20-30/8-12	120/47	61/24	100/39	---
	43/17	40/16	48/19	51-53/20-21	41-53/16-21	17" 20-30/8-12	120/47	63/25	100/39	---
	46/18	40/16	48/19	51-53/20-21	41-53/16-21	18" 20-30/8-12	120/47	66/26	100/39	---
	51/20	40/16	48/19	51-53/20-21	41-53/16-21	20" 20-30/8-12	120/47	71/28	100/39	---

TECHNICAL SPECS

Weight capacity: 136kg / 300lbs
 Max. speed: 7.2kph / 4.5mph
 Range: 40km / 25mi.
 Turn radius: 61cm / 24"
 Ground clearance: 10cm / 4"

Motor: DC 24V, 180W(Rating) / 630W(Max)
 Controller: PG VR2 60A
 Battery: 12V / 22NF x 2pcs
 Charger: 5A off-board
 Gradient: 10°

Caster: 9" (2.80/2.50-4) foam filled tire
 Drive wheel: 14"(3.00-8) foam filled tire
 Brake: Intelligent, regenerative, electromagnetic brakes

Medop cc

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za