

OMEGA A1 POWER



NAME	OMEGAA1 POWER
PRODUCT CODE	A1POW
SIZE	12", 14", 22'
FRAME	MILD STEEL FOLDING
BATTERIES	2X36AH/12 VOLT
MOTOR	250W
CASTORS	8X2" FOAM FILLED GRAY
SPEED & RANGE	6KM/H, 20KM IDEAL CONDITIONS
REAR RIMS & TYRES	121/2 X 21/4 ALUMINIUM RIM 121/2 X 21/4 FOAM FILLED TYRE
ARMS & FOOT	HEIGHT ADJUSTABLE DESK ARM
UPHOLSTERY	BLACK REINFORCED NYLON – TAB
ADJUSTABILITY	HEIGHT ADJUSTABLE BACKREST
GRADIENT	12 DEGREES
CONTROLS	DYNAMIC SHARK – REMD01B POWER MODULE – PMB50
LOAD CAPACITY	LOAD CAPACITY – 120KG GROSS WEIGHT – 65KG
PACKING	BATTERIES – 25KG BOX 1 – 34X80X82CM BOX 2 – 32X35X27CM



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za