

BASIK STANDING

PRODUCT CODE ALK158



Model /**ALK158**
Name /*Electric Wheelchair*



- 1) Spray powder frame with breathable mesh cushion
- 2) 8" big castor, 12" PU rear wheel
- 3) Detachable plate, the armrest can be turned backwards
- 4) can go up and down and walk upright by electric motor
- 5) with safe light
- 6) Seat width of 40cm

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za