

# CTM REHAB POWER (STANDARD, WITH TILT, BED TYPE)



*Lights & cushion not included*



## Features

- \* Maximum Speed: 8kmph/5mph
- \* Weight Capacity:  
7200: 7200T: 162kg / 357lbs  
7200B: 136kg / 300lbs
- \* Suspension System
- \* Standard 75ah Batteries and Easy Accessible Design
- \* Manual Tilting and Adjustable Back Angle
- \* Height Adjustable / Flip Up Armrests

<b>NAME</b>	CTM REHAB POWER
<b>PRODUCT CODE</b>	HS7200T
<b>SEAT WIDTH</b>	16", 18", 20"
<b>BACKREST</b>	ANGLE ADJUSTABLE
<b>CASTOR</b>	9" PNEUMATIC
<b>REAR WHEEL</b>	14" PNEUMATIC
<b>BATTERIES</b>	SEALED 50 A/H
<b>ARMS</b>	FLIP BACK HEIGHT & WIDTH ADJUSTABLE
<b>FOOTREST</b>	REMOVABLE
<b>FRAME</b>	FIXED MILD STEEL
<b>OPTION</b>	TILT IN SPACE

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za