

STORM 4

The *Invacare Storm⁴* powerchair was created to satisfy complex requirements regarding **configurability, adaptability and functionality**.

The *Storm⁴* has impressive outdoor and indoor performance providing a safe and stable drive experience without compromising on manoeuvrability. Thanks to its contemporary design, it also ensures the powerchair is always elegant and stylish. Key features include:

- **Highly configurable**
- **Easy handling**
- **Aesthetic design**
- **LiNX Smart Technology**



Technical data

For more comprehensive clients information about this product, including the product's user manual, please visit www.invacare.co.uk

	Seat width	Seat depth	Backrest height	Total height	Backrest angle	Armrest height	Total weight
<i>Modulite™</i>	380 - 530 mm	410 - 510 mm	480, 540, 560 mm ^{1,10}	1020 - 1190 mm	90 - 120°	275 - 340/ 325 - 390 mm ¹¹ 245 - 310/ 295 - 360 mm ¹²	174 kg ³
<i>Recaro</i>	490 - 530 (360) mm	460 - 510 mm	770 - 830 mm ¹	1300 mm	90 - 135°	250 - 340 / 290 - 380 mm	174 kg ³
	Seat height	Length of footrests	Seat angle	Driving unit width	Total length inc. footrests	Total length without footrests	Max. user weight
<i>Invacare Storm⁴</i>	400 - 650 mm ⁴ 450 / 480 mm without lifter	290 - 460 mm ⁵	0° - 25° ⁹	630 mm (wheelbase)	1190 mm ⁷	910 mm	150 kg
	Battery capacity	Turning diameter	Max. climbable obstacle height	Max. safe slope	Motor power	Speed	Driving range
<i>Invacare Storm⁴</i>	73.5 Ah	1770 mm	60 / 100 mm	20% (11.3°) 15% (8°) with lifter	340 W	6 / 10 km/h	up to 36 km ⁸



Seat tilt

Manual settable or powered tilt with maximum 25° for optimum comfort and posture.



Telescopic seat frame

This frame is the heart of the *Modulite* seat design, supporting a number of positioning options.



Driving unit design

Attractive and modern design with easy to change, colourful shroud covers.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za