


TDX SP2^{NB} POWER

The remarkable *Invacare* **TDX SP2^{NB}**, with its **compact narrow base** boasts **durability** and **ultimate performance**. Great manoeuvrability indoors and out, with a wide choice of seating for **optimum comfort and support** makes the **TDX SP2^{NB}** the powerchair for individuals who need full functionality, but do not want to compromise on style. Key features include:

- Enhanced stability
- Intuitive driving experience
- Simply Smart LiNX control system



Technical data

	Seat width	Seat depth	Seat height	Backrest height	Armrest height	Length of footrests	Seat angle
Modulite™ 1-piece seat plate	440 - 480 mm ¹	410 - 510 mm	450 mm ⁴	480 mm ⁷	275 - 340 mm	290 - 460 mm ⁸	0 - 8°
Modulite™ telescopic seat	380 - 530 mm ²	410 - 510 mm 510 - 560 mm ³	450 / 475 mm ^{4,5} 475 - 775 mm ⁶	480, 540, 560 mm ⁷	245 - 310 mm	290 - 460 mm ⁸	0 - 20° / 0 - 30° ⁵
	Backrest angle	Total height	Driving unit width	Total length inc. Footrests	Total length without footrests	Max. User weight	Battery capacity
Modulite™ 1-piece seat plate	90° - 120°	1010 mm ¹¹	610 mm ¹⁰	1160 mm ⁹	980 mm	150 kg	50 Ah
Modulite™ telescopic seat	90° - 120°	1010 - 1170 mm ¹¹					
	Turning diameter	Max. Climbable obstacle height	Max. Safe slope	Motor capacity	Speed	Driving range	Total weight
	1120 mm	75 mm	9° / 16 %	2 x 340 W	up to 26 km ¹²	6 and 10 km/h	138 kg ¹¹



Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za