

CTM HS115 MICRO SCOOTER



HS-115			
Type	Mobility Scooters	Curb Climbing	35 mm / 1.4"
Model Number	HS-115	Turning Radius	1090 mm / 43"
Overall Length	1010 mm / 39.8"	Supension	N/A
Overall Width	520 mm / 20.5"	Brake	Electro-Mechanical
Overall Height	905 mm / 35.6"	Seat Type	Swivel Padded Foldable
Wheels: Front	195 mm / 8"	Seat Width	425 mm / 17"
Wheels: Rear	195 mm / 8"	Motor Size	250W, 4700 r.p.m.
Weight w/ Batteries	37 kg / 82 lbs	Battery Size	(2) 12V. 12Ah
Max. Speed	6.4 kmph / 4 mph	Battery Weight	9.2 kg / 20 lbs
Weight Capacity	115 kg / 255 lbs	Travel Range	10 km / 6.2 Miles
Ground Clearance	30 mm / 1.2"	Battery Charger	2A Off Board
Grade Climbable	6 degree	Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob
Color	Red		

* Subject to change without notice.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za