

# CTM HS118 MICRO SCOOTER



Type	Mobility Scooters	Curb Climbing	35mm/1.4"
<b>Model Number</b>	HS-118	<b>Turning Radius</b>	1270mm/50"
<b>Overall Length</b>	1000mm/39.4"	<b>Suspension</b>	N/A
<b>Overall Width</b>	510mm/20"	<b>Brake</b>	Electro-Mechanical
<b>Overall Height</b>	890mm/35"	<b>Seat Type</b>	Padded Foldable
<b>Wheels: Front</b>	195mm/8"	<b>Seat Width</b>	425mm/17"
<b>Wheels: Rear</b>	195mm/8"	<b>Motor Size</b>	250W, 4700 rpm
<b>Weight w/ Batteries</b>	40kg/88lbs	<b>Battery Size</b>	(2) 12V • 12Ah
<b>Max. Speed</b>	6.4kmph/4mph	<b>Battery Weight</b>	9.2kg/20lbs
<b>Weight Capacity</b>	115kg/253lbs	<b>Travel Range</b>	10km/6.2miles
<b>Ground Clearance</b>	30mm/1.2"	<b>Battery Charger</b>	1.8A Off Board
<b>Grade Climbable</b>	8 degrees	<b>Electronics</b>	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
<b>Colour</b>	Burgundy, Blue		

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za