

# CTM HS265 SCOOTER



HS-265			
Type	Mobility Scooters	Curb Climbing	50 mm / 2"
Model Number	HS-265	Turning Radius	1060 mm / 41.7"
Overall Length	1020 mm / 40.2"	Supension	N/A
Overall Width	500 mm / 20"	Brake	Electro-Mechanical
Overall Height	890 mm / 35"	Seat Type	Swivel Padded Foldable
Wheels: Front	200 mm / 8"	Seat Width	425 mm / 16.7"
Wheels: Rear	200 mm / 8"	Motor Size	200W, 4650 r.p.m
Weight w/ Batteries	40.1 kg / 88.4 lbs	Battery Size	(2) 12V. 12Ah
Max. Speed	8 kmph / 5 mph	Battery Weight	9.5 kg / 20.9 lbs
Weight Capacity	136 kg / 300 lbs	Travel Range	11 km / 6.8 Miles
Ground Clearance	40 mm / 1.6"	Battery Charger	2A On Board
Grade Climbable	8 degree	Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob
Color	Burgundy - Blue		

\* Subject to change without notice.

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za