

## CTM HS295 AND HS295U SCOOTER



### FEATURES:

Maximum speed: 8kmph  
 Capacity: 136kg  
 Wireless design for trouble free dismantling  
 Combination battery and charger pack  
 Standard Delta handlebar

### TECHNICAL SPECIFICATIONS

Overall Length	1010 mm / 39.8"
Overall Width	550 mm / 21.7"
Overall Height	900 mm / 35.4"
Front Wheels	210 mm / 8 1/4 "
Rear Wheels	210 mm / 8 1/4 "
Weight W/ Batteries	50.3 kg / 110.8 lbs
Max. Speed	8 kmph / 5 mph
Weight Capacity	130 kg / 287 lbs
Ground Clearance	60 mm / 2.4"
Grade Climbable	8 degree
Curb Climbable	40 mm / 1.6"
Turning Radius	1300 mm / 51.2"
Brake	Electro-Mechanical
Seat Type	Swivel Padded Foldable
Seat Width	445 mm / 17.5"
Motor Size	200W, 4650 r.p.m
Battery Size	(2) 12V. 18Ah
Battery Weight	13.4 kg / 29.5 lbs
Travel Range	16 km / 9.9 Miles
Battery Charger	2A Off Board
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob



### FEATURES:

Maximum speed: 8kmph  
 Capacity: 130kg  
 Wireless design for trouble free dismantling  
 Combination battery and charger pack  
 Two 12V.18Ah batteries for longer travel range  
 Standard Delta handlebar

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za