

CTM HS320 MID-LINE SCOOTER



Large Seat Option



Solid Tyre Option

CTM MID-LINE SCOOTER – HS320

Type	Mobility Scooters	Curb Climbing	50mm/2"
Model Number	HS-320	Turning Radius	1020mm/40.2"
Overall Length	1090mm/42.9"	Suspension	N/A
Overall Width	550mm/21.7"	Brake	Electro-Mechanical
Overall Height	870mm/34.3"	Seat Type	Light Weight Foldable Swivel
Wheels: Front	230mm/9"	Seat Width	460mm/18"
Wheels: Rear	230mm/9"	Motor Size	450W 3700 rpm
Weight w/ Batteries	53.8kg/118.7lbs	Battery Size	(2) 12V • 18 Ah
Max. Speed	8kmph/5mph	Battery Weight	12.7kg/28lbs
Weight Capacity	135kg/300lbs	Travel Range	135km/8.4 miles
Ground Clearance	40mm/1.6"	Battery Charger	2A Off Board
Grade Climbable	12 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za