

CTM HS580 SUPREME SCOOTER



CTM SUPREME SCOOTER – HS580

Solid Tyre Option

Type	Mobility Scooters	Curb Climbing	50mm/2"
Model Number	HS-580	Turning Radius	1520mm/59.8"
Overall Length	1280mm/50.4"	Suspension	N/A
Overall Width	610mm/24"	Brake	Electro-Mechanical
Overall Height	960mm/37.8"	Seat Type	Premium Upholstery Mid-Back w/ Headrest Swivel Seat-Sliding Adjustment
Wheels: Front	260mm/10"	Seat Width	455mm/18"
Wheels: Rear	260mm/10"	Motor Size	500W 3400 rpm
Weight w/ Batteries	87.9kg/194lbs	Battery Size	(2) 12V• 36 Ah
Max. Speed	8kmph/5mph	Battery Weight	24.7kg/55lbs
Weight Capacity	135kg/300lbs	Travel Range	30km/18.6 miles
Ground Clearance	40mm/1.6"	Battery Charger	5A Off Board
Grade Climbable	12 degrees	Electronics	On/Off Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za