

CTM HS585 SCOOTER



Solid Tyre Option

HS-585			
Type	Mobility Scooters	Curb Climbing	50 mm / 2"
Model Number	HS-585	Turning Radius	1480 mm / 58.3"
Overall Length	1280 mm / 50.4"	Supension	N/A
Overall Width	610 mm / 24"	Brake	Electro-Mechanical
Overall Height	1110 mm / 43.7"	Seat Type	Swivel Mid-Back W/ Headrest & Seat Sliding Mechanism
Wheels: Front	260 mm / 10"	Seat Width	455 mm / 18"
Wheels: Rear	260 mm / 10"	Motor Size	500W, 3400 r.p.m
Weight w/ Batteries	89.3 kg / 196.9 lbs	Battery Size	(2) 12V. 36Ah
Max. Speed	8 kmph / 5 mph	Battery Weight	24.7 kg / 55 lbs
Weight Capacity	135 kg / 300 lbs	Travel Range	34 km / 21 Miles
Ground Clearance	40 mm / 1.6"	Battery Charger	5A Off Board
Grade Climbable	12 degree	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Color	Burgundy, Silver, Blue, Turquoise		

* Subject to change without notice.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za