

INVACARE LEO

PRODUCT CODE: 1579613



Features and Accessories



Finite tiller movement
Users can easily adjust the tiller to suit their needs.



Brake light
Rear brake lighting warns other people that the scooter is slowing down, thus preventing collisions.



Comfort seating
Swivel functionality and sliding rails ensure optimal seating position.



Two-step disengaging lever
This unique concept prevents the scooter from accidentally freewheeling.

Technical data

470 mm	410 mm	460 - 560 mm	475 mm	590 mm	1220 mm
59 kg (excl. batteries) 83 kg (incl. batteries)	136 kg	2 x 36 Ah	2620 mm	60 mm	70 mm
10°	240 W	38 km	8 km/h		

* Seat in maximum fifth rearward locking position from the front.

Colour of shrouds



Moonstone Silver Onyx Blue Ruby red

Upholstery



Black vinyl

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za